GEBRAUCHSANLEITUNG
Originalversion
OPERATING MANUAL
MODE D'EMPLOI
MANUAL DE INSTRUCCIONES

CT 72/2 | CT 72/4 | CT 72/2 TT

DURCHSICHT - THERMOSTAT
TRANSPARENT THERMOSTAT
THERMOSTAT TRANSPARENTS
TERMOSTATOS TRANSPARENTES

SI Analytics

a **xylem** brand

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Operating Manual..... Page 53 ... 102

Important notes:

The operating manual is part of the product. Before initial operation, please carefully read and observe the operating manual and keep it. For safety reasons the product may only be used for the purposes described in these present operating manual. Please also consider the operating manuals for the devices to be connected.

All specifications in this operating manual are guidance values which are valid at the time of printing. However, for technical or commercial reasons or in the necessity to comply with the statutory stipulations of various countries, the manufacturer may perform additions to the product without changing the described properties. A potentially more recent version of this manual is available on our internet website. The German version is the original version and binding in all specifications!

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TABLE OF CONTENT

QUALITY ASSURANCE	55
GUARANTEE.....	55
EC – DECLARATION OF CONFORMITY	55
YOUR CONTACTS AT XYLEM ANALYTICS GERMANY	55
UNPACKING AND CHECKING	55
INTENDED USE	56
Description	56
OPERATOR RESPONSIBILITY – SAFETY RECOMMENDATIONS	56
OPERATOR RESPONSIBILITY.....	57
SAFETY NOTES.....	58
Explanation of safety notes.....	58
Explanation of other notes.....	58
Safety recommendations.....	59
Recycling and Disposal	61
1. TECHNICAL DATA OF TRANSPARENT THERMOSTATS CT72/2, CT72/4, CT72/2 TT	62
1.1. Technical data of immersion thermostat CT 72.....	62
1.2. Technical data of transparent baths.....	63
2. OPERATING CONTROLS AND FUNCTIONAL ELEMENTS.....	64
2.1. Overall view of transparent thermostats CT 72/2, CT 72/4, CT 72/2 TT.....	66
3. PREPARATIONS	67
3.1. Bath liquids.....	67
3.2. Filling / draining	67
3.3. Modification for using the measuring stand AVS or manual measurement inserts	68
3.4. Countercooling	69
3.5. Magnetic stirrer (option).....	70
4. OPERATING PROCEDURES	71
4.1. Power connection.....	71
4.2. Switching on / Start - Stop	71
5. T SETTING OF TEMPERATURES.....	72
5.1. Using the pre-settings in the T menu	72
5.2. Direct setting of temperatures.....	73
6.  SAFETY INSTALLATIONS, WARNING FUNCTIONS	73
6.1. Low level protection.....	74
6.2. Switch-over from warning to shutdown function.....	75
6.3. Over and Sub temperature warning function	76

7. MENU MENU FUNCTIONS	77
7.1. MENU PROGRAM – START	78
7.2. MENU PROGRAM – creation, administration	80
7.3. MENU CONFIG - Configuration of unit	82
7.3.1. REMOTE	83
7.3.2. AUTOSTART	83
7.3.3. Setting of clock and date	84
7.3.4. RESET - Factory settings	84
7.4. MENU CONTROL – Control characteristics and parameters	85
7.4.1. Dynamik internal	85
7.4.2. Control parameters – XP-, TN-, TV- INTERNAL	86
7.5. MENU SERIAL - BAUDRATE, HANDSHAKE, PARITY	87
7.6. MENU ATC - Absolut Temperature Calibration	88
7.6.1. ATC STATUS - YES / NO	90
7.6.2. CALIBRATION TYPE: 1 -/ 2 -/ 3 POINT	90
7.6.3. Example: 3-point calibration for internal control	91
7.7. MENU LIMITS	93
7.7.1. Limits for internal control	93
8. TROUBLESHOOTING GUIDE / ERROR MESSAGES	94
9. ELECTRICAL CONNECTIONS	96
10. REMOTE CONTROL	97
10.1. Setup for remote control	97
10.2. Communication with a PC or a superordinated data system	97
10.3. List of commands	98
10.4. Status messages	100
10.5. Error messages	100
11. CLEANING / REPAIRING THE UNIT	101

Copyright

© 2022, Xylem Analytics Germany GmbH

Reprinting - even as excerpts - is only allowed with the explicit written authorization.

Germany. Printed in Germany.

Quality assurance

Dear customer,

Xylem Analytics Germany implements a **Quality Management System** certified according to DIN EN ISO 9001.

This guarantees the presence of organizational structures which are necessary to ensure that our products are developed, manufactured and managed according to our customers expectations. Internal and external audits are carried out on a regular basis to ensure that our QMS system is fully functional.

Please inform us if, despite our precautionary measures, you should find any product defects. You can thus help us to avoid such faults in future.

Guarantee

We provide guarantee for the device described for two years from the date of purchase. This guarantee covers manufacturing faults being discovered within the mentioned period of two years. Claim under guarantee covers only the restoration of functionality, not any further claim for damages or financial loss.

Improper handling/use or illegitimate opening of the device results in loss of the guarantee rights. The guarantee does not cover wear parts.

To ascertain the guarantee liability, please return the instrument and proof of purchase together with the date of purchase freight paid or prepaid.

EC – Declaration of Conformity

The corresponding declaration of conformity of the device can be found on our homepage. It will also be made available to you on request.

Your contacts at Xylem Analytics Germany

Please get in contact with us (see backside of this manual) or the authorized agent who supplied you with the unit if you have any further questions.

The following specifications should be given when product enquiries are made:

– **Unit name** printed on the front of the unit and as specified on the name plate.

Unpacking and checking

Unpack the circulator and accessories and check for damages incurred during transit. These should be reported to the responsible carrier, railway, or postal authority, and a request for a damage report should be made. These instructions must be followed fully for us to guarantee our full support of your claim for protecting against loss from concealed damage. The form required for filing such a claim will be provided by the carrier.

Intended use

The transparent thermostats CT 72/2, CT 72/4, CT 72/2 TT are intended to temper liquid media (water, silicone oil or mineral oil, ethanol and polyglycols) in a bath tank. Precise information about the bath media and their operating temperature can be found in Section 3.1 ("Bath liquids ") in this operating manual.

The transparent thermostats are not suitable for direct temperature control of foods, semi-luxury foods and tobacco, or pharmaceutical and medical products. Direct temperature control means unprotected contact of the object with the bath medium (bath fluid).

Description

The main functional elements are the heater, circulation pump, and control electronics.

The transparent thermostats CT72 are operated via the splash-proof keypad. The microprocessor technology allows different values to be set, stored, and displayed on the VFD -DISPLAY.

Temperature- and time-dependent processes can be stored and executed using the integrated programmer.

The adjustable PID cascade temperature control automatically adjusts the heat supply to the thermal requirements of the bath.

Absolute Temperature Calibration (ATC3) provides high temperature stability at all points in the bath. With the 3-point calibration, an offset is adjusted at three temperatures to ensure an accurate temperature pattern at the selected spot in the bath over the entire temperature range.

Electrical connections:

RS232 interface for modern process technology without an additional interface.

Alarm output for external alarm message or control of refrigerating baths.

The excess temperature protection according to IEC 61010-2-010 is a safety device independent from the control circuit. Its safety value can be displayed and adjusted on the VFD -DISPLAY

Operator responsibility – Safety recommendations

The products of SI Analytics® ensure safe operation when installed, operated, and maintained according to common safety regulations. This section explains the potential dangers that may arise when operating the circulator and also specifies the most important safety precautions to preclude these dangers as far as possible.

Operator responsibility

Make sure all persons tasked with operating, installing, and maintaining the unit have read and understand the safety information and operating instructions.

Set up, maintenance, repair

Only qualified personnel is authorized to perform configuration, installation, maintenance and repairs of the thermostat.

Appropriate operation:

Routine operation can also be carried out by untrained personnel who should however be instructed by trained personnel.

The personnel operating the units should be regularly instructed about the dangers involved with their job activities as well as measures to avert these dangers.

When using hazardous materials or materials that could become hazardous, the circulator may be operated only by persons who are absolutely familiar with these materials and the circulator. These persons must be fully aware of possible risks.

Handling

Make sure the product is checked for proper condition regularly (depending on the conditions of use).

Regularly check (at least every 2 years) the proper condition of the mandatory, warning, prohibition and safety labels.

Make sure that the mains power supply has low impedance to avoid any negative effects on instruments being operated on the same mains.

Magnetic radiation may affect other devices with components sensitive to magnetic fields (e.g., monitors). We recommend maintaining a minimum distance of 1 m.

Do not expose the unit to sunlight.

Safety notes

These notes are intended to enhance your own safety consciousness. We have set the highest quality standards for ourselves and this unit during development and production. Every unit meets relevant safety regulations.

The correct unit usage and proper handling is however solely your responsibility.

Explanation of safety notes

In addition to the safety warnings listed, warnings are posted throughout the operating manual. These warnings are designated by an exclamation mark inside an equilateral triangle. "Warning of a dangerous situation (Attention! Please follow the documentation)." The danger is classified using a signal word. Read and follow these important instructions for averting dangers.

Warning

Describes a **possibly** highly dangerous situation. If these instructions are not followed, serious injury and danger to life could result.

Caution

Describes a **possibly** dangerous situation. If this is not avoided, slight or minor injuries could result. A warning of possible property damage may also be contained in the text.

Notice

Describes a **possibly** harmful situation. If this is not avoided, the product or anything in its surroundings can be damaged.

Explanation of other notes

Note!

Draws attention to something special.

Important!

Indicates usage tips and other useful information.

This icon is used in the operating instructions to indicate flashing values or parameters which have to be set or confirmed.

Safety recommendations

Caution:

Handling chemicals, compounds and other

During the use of thermostats, viscosimeters are tempered directly inside the thermostat, e.g. submerged. You can also temper test tubes, Erlenmeyer flasks and similar items.

We do not know which substances are contained within these vessels.

Many substances are:

- inflammable, easily ignited or explosive
- hazardous to health
- environmentally unsafe

The user alone is responsible for the handling of these substances!

Our advice:

- If in doubt, consult a safety specialist.
- Read the product manufacturer's or supplier's "EC – SAFETY DATA SHEET"
- Read relevant regulations concerning dangerous materials.
- Observe relevant guidelines for laboratories in your country.

The following measures were taken for the protection of the operator:

- Protection Class I according to VDE 0106 T1 (IEC536) i.e. protection against electric shocks by grounding all parts which carry the risk of electric contact.

Notice:

Check the safety installations at least twice a year!

- Excess temperature protection according to IEC 61010-2-010.
With a screwdriver turn back the adjustable excess temperature protection until the shut-down point (actual temperature).
- Low level protection according to IEC 61010-2-010.
To check the function of the float, it can be manually lowered with a screwdriver for example.

The following measures were taken for the protection of the operator:

- Protection IP 20 according to EN 60529, i.e. regarding the protection against accidentally touching live parts and damage by foreign matter, it has been ensured that foreign bodies with a thickness or diameter of more than 12 mm cannot penetrate.

Notice:

No special precautions were taken against the penetration of water and dust.

- The device should therefore not be used in a dusty atmosphere or in the neighborhood of spray water.
- Do not insert wires or tools in any of the openings.

Notice:

Burn hazard

- Some parts of the device may become extremely hot during continuous operation.
- Always protect yourself against accidental contact.

Before initial operation

These operating instructions must be carefully studied! It contains important information on the connection to the local mains supply, correct unit usage and safe handling.

Check for transportation damage during unpacking.

Do not try to start up a damaged unit before the damage has been repaired or you have ascertained the effect of the damage. Get in contact with supplier and/or carrier for settlement of damage claims

Repairs, alterations or modifications must only be carried out by specialist personnel. Considerable damage can be caused by improper repairs. Our service is at your disposal for repair work.

Never remove warning labels!

Initial operation

Only connect the unit to a power socket with an earthing contact (PE – protective earth). Electrical connections and any other work must be performed by qualified personnel only. The power supply plug serves as a safe disconnecting device from the line and must always be easily accessible. Never operate units with damaged mains power cables.

Place the unit on an even surface on a base made of nonflammable material.

Set the excess temperature safety installation at least 25 °C below the fire point of the bath fluid.

Operation

Only use this unit solely for the intended application.

Never operate the unit without bath fluid in the bath! Only use the heat transfer liquids recommended by us.

Do not operate the unit with wet or oily hands.

Do not expose the unit to spray water or immerse it in water.

Do not drain the bath fluid while it is hot.

Cleaning

Do not clean the unit with solvents (fire risk!) – a wet cloth soaked in household detergent is normally sufficient.

Transportation

Do not move the instrument while it is running or hot. Always empty the bath before moving the unit. Transport the unit with care. Sudden jolts or drops may cause damage in the interior of the unit.

Notice:

Shock hazard

Complete separation from mains is required when:

- all dangers caused by this device are to be avoided,
- cleaning is carried out;
- repairs or maintenance work is about to be carried out

Complete separation means:

Pull out the mains plug!

Recycling and Disposal

Please observe the applicable local or national regulations concerning the disposal of "waste electrical and electronic equipment"

In Germany, disposal of the device in household waste or via municipal collection points is not permitted. **Devices that are no longer used must be sent to the respective manufacturer for disposal with paid postage and the note "FOR DISPOSAL"!** The devices will then be disposed of at the manufacturer's expense.

Detailed information on "Recycling and disposal" can be found on our homepage in the customer support section under "WEEE disposal concept". If you have any further questions about disposal, please contact the manufacturer directly (see the back of these operating manual).

The device and its packaging have been manufactured as far as possible from materials that can be disposed of in an environmentally friendly manner and recycled properly.

i This device contains batteries. Batteries must not be disposed of with household waste. The end user is obliged to dispose of used batteries (including non-hazardous batteries) at a collection point set up for this purpose or at the point of sale for proper recycling.

1. Technical data of Transparent Thermostats CT72/2, CT72/4, CT72/2 TT

Translation of the legally binding German version
(Version as of 18. January 2012)

Country of origin: Germany / Made in Germany

1.1. Technical data of immersion thermostat CT 72

Temperature ranges:

Operating temperature range			CT72/2, CT72/4	CT72/2TT
with additional cooling	°C		+ 5 ... + 150	
with refrigerated circulator (option)				- 40 ... + 150
Working temperature range				
(20 °C ambient temperature)	°C		+ 30 ... + 150	+ 30 ... + 150
with tap water cooling*	°C		+ 20 ... + 150	+ 20 ... + 150
with flow-through cooler CK300/CK310 (option)	°C		+ 5 ... + 150	+ 5 ... + 150
Temperature stability	K		± 0,02	
Temperature drift at changing ambient temperature between 5 and 40 °C:			max. 0,002 K/1°C	
Temperature control			PID3 cascade temperature control	
<u>Heater capacity:</u>				
	at 230 V	W	2000	
	at 115 V	W	1000	
<u>Pump</u> :				
	Pressure max. at 0 liter	bar	0.35	
	Flow rate at 0 bar	l/min	15	
<u>Connections:</u>				
electrical:				
connector: control cable to the flow-through cooler CK300/CK310		Vdc/mA	24-0 / max. 25	
Computer interface			RS232	
mechanical				
nozzles for flow-through cooler. CK300 / CK310			2 x Ø 12 mm	
nozzles for tap water cooling			2 x Ø 8 mm	
Overall dimensions	(WxDxH)	cm	13x16x33	
Weight		kg	5.0	
Mains power connection	230 V/50-60 Hz	V/ Hz	190 ... 253 / 50-60	
Current draw	(at 230 V)	A	9	
Mains power connection	115 V/60 Hz	V/ Hz	90 ... 127 / 60	
Current draw	(at 115 V)	A	9	

* If tap water is used, an operating temperature of approx. 3 °C above the given tap water temperature can be maintained.

Safety installations according to IEC 61010-2-010:

Excess temperature protection	adjustable from 0 °C ... 230 °C
Low liquid level protection	float switch
Classification according to DIN 12876-1	class III

Supplementary safety installations

Early warning system for low level	float switch
High temperature warning function	optical + audible (in intervals)
Low temperature warning function	optical + audible (in intervals)
Supervision of working sensor	plausibility control
Reciprocal sensor monitoring between working and safety sensors	difference >35 °C
Alarm message	optical + audible (permanent)
Warning message	optical + audible (in intervals)

Environmental conditions according to IEC 61 010-1:

Use only indoor.

Altitude up to 2000 m - normal zero.

Ambient temperature: +5 ... +40 °C

Air humidity:

Max. rel. humidity 80 % for temperatures up to +31 °C,

linear decrease down to 50 % relative humidity at a temperature of +40 °C

Max. mains fluctuations of ±10 % are permissible.

Protection class according to IEC 60 529	IP20
The unit corresponds to Class I	
Overvoltage category	II
Pollution degree	2

Caution:

The unit is not for use in explosive environment

Standards for interference resistance according to EN 61326-1

This unit is an ISM device classified in Group 1 (using high frequency for internal purposes) Class A (industrial and commercial range).



1.2. Technical data of transparent baths

		CT72/2, CT 72/2 TT	CT72/4
<u>Filling volume:</u>	L	ca. 15	ca. 27
<u>Overall dimensions</u>	(WxDxH) cm	35,5 x 25 x 37	60,5 x 25 x 37
<u>Weight:</u>	(without bath liquid) kg	Approx. 13,5	Approx. 28,0
<u>Material:</u>		High-grade steel, Insulating glass (mineral glass)	
<u>Bath liquid:</u>		Water, Silicone oil, Mineral oil, Polyglycole	
<u>Country of origin:</u>		Made in Germany	
<u>Temperature range:</u>			
CT72/2, CT72/4	+5 ... +150°C	With built-in cooling coil in combination with a flow-through cooler CK300/CK310 (option) or a refrigerated circulator (option)	
CT72/TT	-40 ... +150°C	With special built-in cooling coil TT in combination with a refrigerated circulator (option)	

2. Operating controls and functional elements

Front view CT 72

Rear view CT 72

1

Mains power switch, illuminated

2

VFD-DISPLAY

Header: Control indicators see sections 2.1 and 2.2

Line 1: Actual value internal

Line 2: Working temp. setpoint, constantly S xxx.xx

Line 3: Actual value (I = internal)
Alternating with the display in line 1

2.1

Control indicators in the header:

Heating / Cooling / Alarm /
Remote control

2.2

Control indicators in the header:

Temperature indication **Internal** actual value.
(**External** not possible on this unit)

Temperature indication in **°C**
(**°F** not possible on this unit)

3

Adjustable excess temperature protection according to IEC 61010-2-010

4		Navigation keys
4.1		<p>1. Key: > ENTER < Start / Stop (pump / heater)</p> <p>2. > ENTER < in the menu Menu item / select submenu for setting Save set value Save selected parameter</p> <p>A beep signals the end of setting</p>
		<p>After the actions Start, Stop and change from VFD Display to standard display the key  is locked for a short time. The above graph “front side” shows an example for standard display.</p>
4.2		<p>1. Key: >Return< Stop (pump / heater)</p> <p>2. >Return< in the menu one menu level down Correction function for parameters or values (prior to ENTER)</p>
		immediately back to standard display
		 -  icon for „keep key pressed down“.
4.3	 	<p>1. Key: >Up / Down <temperature – increase/decrease setpoint Push key quickly for single steps, Keep key pressed for fast change.</p> <p>2. >Up/Down< in the menu selection of menu items / parameters</p>
		Menu keys
4.4		Key: start the menu > warning and safety values<
4.5		Key: start the menu >temperature setpoints<
4.6		Key: display of MENU structure
5	 RS232	Interface RS232: remote control via personal computer
6		5-pin connecting socket for flow-through cooler CK300/CK310 via adapter cable
7		Mains power cable with plug
8		Mains circuit breakers (resettable) 15 A

2.1. Overall view of transparent thermostats CT 72/2, CT 72/4, CT 72/2 TT


The Transparent Thermostats are supplied ready for operation and consist of the following components depending on the ordered model:

- CT 72/2 consisting of: immersion thermostat CT 72 and transparent bath 2
- or CT 72/2 TT consisting of: immersion thermostat CT 72 and transparent bath 2 TT
- or CT 72/4 consisting of: immersion thermostat CT 72 and transparent bath 4

CT 72/2

CT 72/2 TT

CT 72/4

3. Preparations

3.1. Bath liquids

Caution

Carefully read the safety data sheet of the bath liquid used, particularly with regard to the fire point!

If a bath fluid with a fire point of ≤ 65 °C is used, only supervised operation is possible.

Water: The quality of water depends on local conditions.

- Due to the high concentration of lime, hard water is not suitable for temperature control because it leads to calcification in the bath.
- Ferrous water can cause corrosion - even on stainless steel.
- Chloric water can cause pitting corrosion.
- Distilled and deionized water is unsuitable. Their special properties cause corrosion in the bath, even in stainless steel.

Xylem Analytics Germany recommends:

1. Water, in the mixing ratio: 70 % deionized water and 30 % tap water
Temperature range: 5 °C to 80 °C

Important:

Using the Transparent Thermostat CT 72/2, CT 72/2 TT in the **Viscosity Automatic Sampler AVS Pro** with water as the bath liquid, the recommended temperature range is from 10 °C to 40 °C.



2. The use of silicone oil or mineral oil with a maximum viscosity of 30 mm²/s,
Temperature range: -40 °C to 150 °C
3. Ethanol and Polyglycole
Temperature range: -40 °C bis +10 C°

3.2. Filling / draining

Notice

- Prevent humidity from entering into the immersion thermostat CT72.
- Pay attention to the thermal expansion of bath oil during heating to avoid overflowing of the liquid.
- Always turn off the unit and disconnect the mains cable from the power source before cleaning the unit, or before moving the unit.
- Do not drain the bath fluid while it is hot!
Check the temperature of the bath fluid prior to draining (by switching the unit on for a short moment, for example).
- Always turn off the unit and disconnect the mains cable from the power source before cleaning the unit, or before moving the unit.
- Store and dispose the used bath fluid according to the laws for environmental protection.


Filling

- Recommended maximum filling level with water as the bath liquid is the top edge of the window.
- Recommended maximum filling level with oil as the bath liquid is 20 mm below the top edge of the window.

Draining

- Press the mains switch to turn the transparent thermostat off.
- Disconnect the power plug from the mains socket!
- Slide a short piece of tube onto the drain port (17) and hold it into a pail.
- Unscrew the drain tap (18) and empty the unit completely.
- Tighten the drain tap.

3.3. Modification for using the measuring stand AVS or manual measurement inserts

(Contents of delivery)

Bath lid B suitable for 2 measuring stands AVS

or

2 manual measurement inserts for capillary viscometers with holders

CT 72/2, CT 72 TT

- 1 Lay centering insert for measuring stand AVS on the bottom of the bath and push it in direction to edges **A**.
- 2 Lay the clamp on the rear panel and tighten both screws to fix the centering insert to the bottom of the bath.
- 3 The insert for the bath lid is attached with six screws.
- 4 Each of the two openings can accept one measuring stand AVS
- 5 If necessary, immerse manual measurement inserts for visco-meters in the opening.
Immerse viscometers with holders as per the instructions in the corresponding operating manual for capillary viscometers.
- 6 Cover the unused openings with the lid.

CT 72/4:

Bath lid B suitable for 2 measuring stands AVS

and

2 manual measurement inserts for capillary viscometers with holders


3.4. Countercooling

Notice

- Observe the laws and regulations of the water distribution company valid in the location where the unit is operated.
- Make sure that the tubes are securely attached.
- Employ suitable connecting tubing (temperature range).
- Avoid mechanical work (torsion, bending) to be applied to tubing.
- Regularly check the tubing for material defects (e.g. for cracks).

For applications near the ambient temperature, the cooling coil must be connected to the water mains.

Connection to cooling (tap) water

- Using tubing, connect the cooling coil to the tap water supply, and lead the tap water in a sink through the return connector.
- A specific water flow rate of 100 ml/minute is sufficient to compensate for the characteristic temperature.
Recommendation: Use a flow-rate regulator.

15 nozzles 2 x 8 mm dia. for tap water cooling

External Cooling Devices

Recommended in combination with CT72/2, CT72/4

With the SI Analytics® flow-through cooler CK300/CK310, the heat transfer liquid can be cooled down and the transparent thermostat can be rendered independent of tap water.

- Installation and use are described in the instruction manual for the cooler.

16a Pump connector: Feed - OUT

16b Pump connector: Return - IN
nozzles 2 x 12 mm dia. for flow-through cooler (CK300/CK310)


14 5-pin socket for control cable to the flow-through cooler (CK300/CK310)

Recommended in combination with CT72 TT

Using an appropriate refrigerated circulator, the bath liquid can be cooled down to -40 °C.

- Installation and use are described in the instruction manual for the refrigerated circulator.

15 Using the nozzles 2 x 8 mm dia. connect the tubing onto the pump nozzles of the refrigerated circulator.

Insulation of the bath cover B

The insulation reduces condensation of moisture and prevents the formation of ice on the bath cover.

- ❶ Insert a measuring stand AVS in each of the two openings and place it on the bath bottom
or
Immerse manual measurement inserts for viscometers in the opening and insert viscometers with holders.
Cover the unused openings with the lid.
- ❷ Lay on insulation for the bath cover.
- ❸ Lay the metal bath cover on the insulation.

3.5. Magnetic stirrer (option)

Transparent baths CT 72/2, CT 72 TT, CT 72/4:

When using AVS-measuring stands, each measuring point can be equipped with a magnetic stirrer later.

On bottom of the bath there are removable covers with leading-in tubes for the necessary supply cables. These covers are accessible from below.

4. Operating procedures

4.1. Power connection

Caution

- Only connect the unit to a power socket with earthing contact (PE – protective earth)!
- The power supply plug serves as safe disconnecting device from the line and must be always easily accessible.
- Never operate equipment with damaged mains power cables.
- Regularly check the mains power cables for material defects (e.g. for cracks).
- We disclaim all liability for damage caused by incorrect line voltages!

Check to make sure that the line voltage matches the supply voltage specified on the identification plate. Deviations of $\pm 10\%$ are permissible.

4.2. Switching on / Start - Stop

Switching on:

- Turn on the mains power switch (1).

i The unit performs a self-test. Then the software version appears. The display „OFF“ or „R OFF“ indicates the unit is ready to operate.

i The circulator enters the operating mode activated before switching the circulator off:

keypad control mode (manual operation)

or

remote control mode (operation via personal computer).

Start:

- Press  key.
The actual bath temperature is displayed on the VFD COMFORT-DISPLAY.
The circulating pump starts with a slight delay.

Stop:

- Press  key.
- or
- Keep  key pressed.

The unit also enters the safe operating state "OFF" or "r OFF" after a mains power interruption. The temperature values entered via the keypad remain in memory. With the circulator in keypad control mode, press the start/stop key to restart operation. With the circulator in remote control mode, the personal computer must first resend the parameters set via the interface before the circulator may be restarted.

5. **T** Setting of temperatures

5.1. Using the pre-settings in the **T** menu

Press the **T** key to call up the menu for temperature selection.

3 different working temperatures can be adjusted. Their values are freely selectable within the operating temperature range.

- ① The temperatures can be set in start or stop mode.
- ① Press **ESC** key if a value is to be retained

Setting of working temperature in the **T** menu

1. Press the key **T**. The value flashes ✖
2. Select SETPOINT 1 or 2 or 3 using the key ▼ or ▲.
3. Confirm by pressing the ↵ key.

Factory setting:

SETPNT 1 25 °C
SETPNT 2 37 °C
SETPNT 3 70 °C

- ① The circulator uses the new working temperature value for temperature control.

Example: Adjustment/modification of the pre-setting of "SETPOINT 3"

1. Press the **T** key.
2. Select SETPOINT 3 by pressing the ▼ key.
Example: SETPNT 3 / 70.00 °C
3. Keep the ↵ key pressed until the integer digits flash ✖.
(example: <70>)
4. Adjust value by pressing the ▼ key and the ▲ key to 85.00 °C and confirm by pressing the ↵ key.
The decimal digits flash ✖ and can be adjusted if desired.
Confirm once more by pressing the ↵ key.
Example on the left: SETPNT 3 / 85.00.

88.888
SETPNT
3 70.00 ✖

88.888
SETPNT
3 85.00 ✖


SSS
70.93
5 85.00

- ① If the active setpoint (SETPNT) is changed, the new value is immediately used for the control of the working temperature.
The heater control indicator flashes.
- ① If the other two setpoints (not activated for control) are changed the **T** MENU has to be left by pressing the **ESC** key after the decimal digits have been confirmed

Notice: Refer to SETPOINT MAX / MIN in chapter 7.7. MENU LIMITS

5.2. Direct setting of temperatures

The circulator uses the setpoint of SETPNT 1 or 2 or 3 for temperature control

The indicated setpoint temperature can be changed directly any time.

Example: change 25.00 °C to 45.00 °C

1. By pressing the key the circulator switches to the active SETPOINT< example on the left: >SETPNT / 1 25.00°C<. The integer digits flash (example: <25>).

2. By pressing the keys and the value is changed to 45.00 °C and is confirmed by pressing the key. The decimal digits flash and can be adjusted if desired. Confirm once more by pressing the key.

- ① The circulator uses the new working temperature value for temperature control.
- ① The temperatures can be set in start or stop mode.

6. Safety installations, warning functions

Check the safety installations at least twice a year!

Refer to section “Safety recommendations”

SECVAL
(Security Values)

- SAFETMP
- AL-TYPE
- OVERTMP
- SUBTEMP

Settings for the excess temperature protection > **SAFETMP**<

and for the warning functions for high > **OVERTMP**< and low > **SUBTEMP**< temperature are made in a menu which is called up by pressing the key .

Menu item > **AL-TYPE**< allows choosing between a warning and

an alarm cut-off for the menu items > **OVERTMP**< and > **SUBTEMP**<.

6.1. Excess temperature protection

Warning

The excess temperature protection **must** be set at least 25 °C below the fire point of the bath fluid used!

In case of wrong setting there is a fire hazard!

We disclaim all liability for damage caused by wrong settings!

This excess temperature protection is independent of the control circuit. When activated heater and circulating pump are completely shut down.

The alarm is indicated by optical and audible signals (continuous tone) and the error message "ALARM-CODE 14" appears on the VFD COMFORT-DISPLAY together with the ticker:

> EXCESS TEMPERATURE PROTECTOR ALARM-CHECK ADJUSTMENT <

Setting range: 20 °C ... 230 °C

① Rough setting can be effected by using the temperature scale.


Exact setting:

1. Press the key  to display menu >SAFETMP<.

2. Press the  key and the set shutdown value is indicated.

Set the new shutdown value within 30 seconds using a screwdriver. The value is indicated on the VFD COMFORT-DISPLAY

Example: SAFETMP / 50 °C

Recommendation:

Set the excess temperature protection at 5 °C to 10 °C above the working temperature setpoint.

6.1. Low level protection

This low level protection is independent of the control circuit.

If the low liquid level protection device is triggered, a complete shutdown of the heater and circulating pump is effected.

A continuous alarm sounds and together with the ticker: > LOW LEVEL ALARM-FILL MEDIUM < a message appears on the VFD -DISPLAY: "ALARM-CODE 1"

Turn off the unit with the mains switch, refill bath fluid and turn the unit on again!

Notice

When adding bath fluid, always use the type of fluid which is identical with the fluid in the bath.

6.2. Switch-over from warning to shutdown function

SECV AL
AL - TYPE

Factory setting:
>WARNING<

AL - TYPE
WARNING ✱


AL - TYPE
ALARM

If a shutdown of functional elements (e.g. heater, circulating pump) is required when the limit values are exceeded or undercut the circulator can be changed over from warning function >WARNING< to shutdown function >ALARM<.

1. Press the key .
2. Select the menu >SECV AL -AL-TYPE< by pressing the  key.
3. Press the  key and the set parameter will flash ✱✱.
(Example: WARNING)
4. Change the parameter by pressing the  key and confirm by pressing the  key.
or
press the **ESC** key if the parameter is to retained.

ⓘ Setting >WARNING<

A mere warning function with optical and audible warning signal (interval tone) A message appears on the VFD COMFORT-DISPLAY:

	or	
--	----	--

• Setting >ALARM<

Temperature limit with shutdown of heater and circulating pump.
An audible alarm sounds (continuous tone) and a message appears on the VFD COMFORT-DISPLAY:

	or	
---	----	---

6.3. Over and Sub temperature warning function


Over temperature

OVERTMP
200.00

Sub temperature

SUBTEMP
-99.00

If the observance of a working temperature value >SETP< has to be supervised for a sensitive temperature application, then set over and sub temperature warning values. In the example below the SETPOINT 45 °C is surrounded by the values OVERTMP 48 °C and SUBTEMP 43 °C. The electronics immediately register if the actual temperature breaches one of the set limit values. The resulting reaction is defined in a further menu item. (See chapter 6.2. Switch-over from warning to shutdown function)

1. Press the key .
2. By pressing the  or  key select the menu >OVERTMP< or >SUBTEMP<.
3. Press the  key. The integer digits flash
4. Change the values to 87. °C and/or 83. °C by pressing the  and  key and confirm with the  key.
The decimal digits flash and can be adjusted if desired.
Confirm once more by pressing the  key.
See above examples.

 The warning functions are only activated if the actual bath temperature remains within the set limit values for 3 seconds after switch-on.

Recommendation:

Set the over temperature warning value >OVERTMP< 5 °C to 10 °C above the working temperature setpoint.

Set the sub temperature warning value >SUBTEMP< 5 °C to 10 °C below the working temperature setpoint.

7. MENU Menu functions

The term „Menu functions“ refers to settings such as

1. Open the menu by pressing the **MENU** key.
2. Use the   keys to scroll in menu level 1.
3. Press the  key to change to menu level 2.
Press the **ESC** key if settings are to be retained.

MENU ↓

Menu level 1

MENU
P-START 
 

Start program

P-START
STEP

Page 78

MENU
PROGRAM

Administration and creation of programs

Page 80

MENU
CONFIG

Configuration of the unit

- REMOTE – on / off (remote control via RS232)
- AUTOST – AUTOSTART on / off
- TIME / DATE – setting time and date
- RESET – factory settings

Page 82

MENU
CONTROL

Control characteristics and parameters

C-TYPE – Internal or external control

DYNAMICS - internal

Control parameter - XP-, TN-, TV- INTERNAL

Page 85

MENU
SERIAL

Adjustable interface parameters

BAUDRATE, PARITY, HANDSHAKE

Page 87

MENU
ATC

ATC - Absolute Temperature Calibration,

Sensor calibration INTERNAL SENSOR,
3-point calibration

Page 88

MENU
LIMITS

Limitations of temperature and capacity

SETPOINT MAX / MIN - Maximum and minimum setpoint

HEAT MAX – Set maximum heating

COOLING MAX – Set maximum cooling

INTERN MAX / MIN – Limitation of the temperature range

BAND HIGH / LOW – Band limit

Page 93

7.1. MENU PROGRAM – START

Start-Menu

This menu will start a previously set program.

i Requirements:

1. Create a program. (refer to next chapter)
2. Return to the Start-MENU and confirm the desired setting of each MENU item with the key
3. Set a start time (>TIME< >DATE< >YEAR<) if the program is to be started by the internal timer.

Menu level 1


- > STEP< Program start at section 1 ... 10
- > RUNS < Number of repetitions 1 ... 99
- > END< Status at end of program (STDBY/SETPNT)
Standby or last setpoint
- > GO < Time of start (NOW/TIMER)


Level 2

Parameter level


i Press the **ESC** key if a parameter is to be retained.
Correction function for parameters or values (prior to Enter)


- Set program step with and
- example: STEP 1

- Set number of runs with and
- example: 1 run








- Set desired parameters with and .
- i** (STDBY / SETPNT)
Parameter **STanDBY**: the circulator switches to –OFF-.
Parameter **SETPoiNT**: the circulator constantly keeps the temperature at the value of the last step.


- Confirm >NOW< with the key and the program will start immediately
- i** or
start at the set time under parameter (TIMER).
Set time in the example below:
09. August 2009, 11:15 hrs

 <p>Submenu <u>TIMER</u></p>	<p>Parameter level</p>	<p>ⓘ set the time for the start of the program in the submenu >TIMER<.</p>
		<p>>TIME< hours/minutes (hh:mm), set both values one after the other and confirm</p> <ul style="list-style-type: none"> hours flash, set by pressing   +  minutes flash, set by pressing   + 
		<p>>DATE< day/months (TT/MM), set both values one after the other and confirm.</p> <ul style="list-style-type: none"> day flashes, set by pressing   +  month flashes, set by pressing   + 
		<p>>YEAR< year</p> <ul style="list-style-type: none"> Set the year with   and .
		<p>ⓘ The program starts at the set time.</p>
	<p>ⓘ Display of time until start: In line 3 the notice >TIMER< and the set values for „TIME“ and „DATE/YEAR“ are alternately indicated</p> <p>ⓘ Check the correct setting of the internal real time clock if required (see MENU CONFIG)</p>	

The started program

After the start the program will indicate the currently calculated setpoint in line 2 S XX.XX. The value increases within the time period >TSLICE< until the target temperature >SETPNT< of the section is reached. If the time period in a section is set to „0“, the next section will not begin until the target temperature has been reached.

Use the edit keys   to scroll to line 3. The display changes approximately every 4 seconds between the current section (STEP XX) and the

- A remaining time of the section
- B remaining time of the program
- C current bath temperature
- I xxx.xx - internal actual value or
- E xxx.xx – external actual value

C2	S XXXX E XXXXX	D RUN – the program has started or PAUSE – the progress of the program has been interrupted by pressing the ESC key. While the time is stopped the temperature will constantly remain at the last calculated setpoint Continue with the ESC key.
D1	S XXXX RUN	
D2	S XXXX PAUSE	

Termination / Interruption of a program

- ❶ The program can be terminated any time by pressing the key **ESC** .
- ❷ In case of power failure the program is interrupted. The circulator switches to **-OFF-**.
- ❸ If the AUTOSTART-function is activated the programmer starts again at the point in time approx. 5 minutes prior to the interruption. However, an uncontrolled change of the bath temperature has occurred.

7.2. MENU PROGRAM – creation, administration


Menu level 1

MENU
PROGRAM

1 program
10 sections

The integrated programmer permits fast and easy programming of setpoint temperature sequences. This temperature sequence is called program. A program is composed of individual sections (STEP). The sections are defined by duration (TSLICE) and target temperature. The target temperature is the setpoint (SETPNT), which is achieved at the end of a section. The programmer calculates the temperature ramp from the difference in time and temperature.

STEP (Nr.)	1	2	4	6	8
SETPNT (°C)	100	180	180	75	75
TIME (hh.mm)	00:20	00:10	00:20	00:50	00:20








- ❶ Sections without set value and time are skipped. They can be defined retroactively and the integrated into the program.

Menu functions

<p>Menu level 1</p> <p style="text-align: center;">▼ ▲</p> <p style="text-align: center;">MENU PROGRAM</p> <p style="text-align: center;">← →</p>	<p style="text-align: center;">>EDIT< Create, administer program</p> <p style="text-align: center;">>STEP< Program step (1 ... 10)</p> <p style="text-align: center;">>SETPNT < Temperature setpoint of step ...</p> <p style="text-align: center;">>TSLICE< Duration of step ...</p> <p style="text-align: center;">>DELETE< delete program step (01 ... 10, ALL)</p> <p>ⓘ Press ESC key, if a parameter is to be retained. Correction function for parameters or values (prior to ENTER)</p>	
<p>Level 2</p> <p style="text-align: center;">PROGRAM EDIT</p> <p style="text-align: center;">← →</p> <p style="text-align: center;">▼</p>	<p>Level 3</p> <p style="text-align: center;">EDIT STEP</p> <p style="text-align: center;">← →</p> <p style="text-align: center;">EDIT 01 STEP</p> <p style="text-align: center;">←</p> <p style="text-align: center;">▼ ▲</p>	<p>Parameter level</p> <p style="text-align: center;">STEP</p> <p style="text-align: center;">1 ✖ (STEP 1 ... 10)</p> <ul style="list-style-type: none"> Set program step with ▼ ▲ and ← ← (Example: EDIT STEP 01) ⓘ For STEP 01 the values for SETPOINT 01 and TSLICE 01 are set one after the other
	<p style="text-align: center;">EDIT 01 SETPNT</p> <p style="text-align: center;">← →</p> <p style="text-align: center;">▼ ▲</p>	<p style="text-align: center;">SETPNT</p> <p style="text-align: center;">10000 ✖ (values within working temp. range)</p> <ul style="list-style-type: none"> Integer digits flash, set by pressing ▼ ▲ + ← Decimal digits flash, set by pressing ▼ ▲ + ←
	<p style="text-align: center;">EDIT 01 TSLICE</p> <p style="text-align: center;">← →</p> <p>(time slice)</p>	<p style="text-align: center;">TSLICE</p> <p style="text-align: center;">00.10 ✖</p> <ul style="list-style-type: none"> Set duration by pressing ▼ ▲ and ←
<p style="text-align: center;">▲</p> <p style="text-align: center;">PROGRAM DELETE</p> <p style="text-align: center;">← →</p>	<p style="text-align: center;">DELETE STEP -- ✖</p>	<ul style="list-style-type: none"> ⓘ Delete program Program steps can be deleted individually or entirely. (STEP 01, 02,... 10, ALL). Set parameters by pressing ▼ ▲ and ←

7.3. MENU CONFIG - Configuration of unit

<p>Menu level 1</p>  <p>MENU CONFIG</p>	 <p>① A RESET can be effected only in the >OFF< mode. Switch off circulator by pressing the  key and call up the menu CONFIGURATION.</p>	
<p>Level 2</p>	<p>Parameter level</p>	<p>① Press the key ESC if a parameter is to be retained. Correction function for parameters and values (prior to ENTER).</p>
 <p>CONFIG REMOTE</p>	 <p>REMOTE OFF ✖</p> <p>or</p> <p>REMOTE ON</p>	<ul style="list-style-type: none"> • Switch on and off remote control by pressing   and  ① Control display in the topline  for Remote ① For remote control refer to page 97 ①  Connect RS232 with PC.
 <p>CONFIG AUTOST</p>	 <p>AUTOST OFF ✖</p> <p>or</p> <p>AUTOST ON </p>	<ul style="list-style-type: none"> • Switch on and off Autostart by pressing   and  <p>AUTOSTART on = on AUTOSTART off = off See WARNING page 83</p>
 <p>CONFIG TIME/DT</p>	<p>Level 3</p>  <p>TIME/DT TIME</p>	<p>Parameter level</p>  <p>TIME 11.15 ✖</p> <ul style="list-style-type: none"> • Hours flash, set by pressing   +  • Minutes flash, set by pressing   + 
 <p>TIME/DT DATE</p>	 <p>DATE 08.09 ✖</p>	 <p>DATE 08.09 ✖</p> <ul style="list-style-type: none"> • Day flashes, set by pressing   +  • Month flashes, set by pressing   + 
 <p>TIME/DT YEAR</p>	 <p>YEAR 09 ✖</p>	 <p>YEAR 09 ✖</p> <ul style="list-style-type: none"> • Year flashes, set by pressing   + 

	Level 3	Parameter level
	 	<ul style="list-style-type: none"> Return to factory settings by pressing  <p>RESET returns all set values to the factory setting except for date and time.</p> <ul style="list-style-type: none"> i A RESET can be effected only in the –OFF- mode. i During the message –RUN- all parameters are reset to factory settings.

7.3.1. REMOTE

Factory setting: OFF

The control electronics offer two ways of adjusting a setpoint.

- Adjustment of setpoint using the keypad or the integrated programmer.
- Adjustment of setpoint via the serial interface RS232 using a PC or a superordinated process control system.

i The topline of the VFD-DISPLAY shows a bright „R“ for remote control; - remote control discontinued.

RS232

IMPORTANT: additional measures for remote control

- i** Connect the circulator to the PC using an interface cable.
- i** Check the interface parameters of both interfaces (Thermostat and PC) and make sure they match. (refer to 10.1 «Setup for remote control»)

7.3.2. AUTOSTART

Warning

For supervised or unsupervised operation with the “AUTOSTART“ function avoid any hazardous situation to persons or property

Take care to fully observe the safety and warning functions of the circulator.

Factory settings: OFF

Notice:

The circulator has been configured and delivered by Xylem Analytics Germany in accordance with the NAMUR recommendations. This means for the start mode that the unit must enter a safe operating status after a power failure. This safe operating status is indicated by the message „OFF“ or „R OFF“ on the VFD COMFORT-DISPLAY.

A complete, all-pole shutdown of the main functional elements such as heater and pump motor is effected.

The values set on the circulator remain saved and the unit is restarted by pressing the start/stop key in manual control.

In remote control mode the values need to be resent by the PC via the interface. If such a safety standard is not required, the NAMUR recommendations can be bypassed with the AUTOSTART function thus allowing a direct start of the circulator by pressing the mains switch or using a timer.

7.3.3. Setting of clock and date

The internal real time clock allows starting a program any time. The clock is set to the local mean time (MEZ) at the factory.

- ① If the unit is operated in a different time zone, the clock can be adjusted in this menu.
- ① Change summer/winter time in this menu

7.3.4. RESET - Factory settings

A Reset will return all values to factory setting except for date and time.



- ① A RESET can be effected in the >OFF< mode only.
Switch off the circulator by pressing the key  and call up the menu CONFIGURATION.

7.4. MENU CONTROL – Control characteristics and parameters

<p>Menu level 1</p> <p>▼ ▲</p> <p>MENU CONTROL</p> <p>← →</p>	<p>The Thermostat is qualified for internal temperature control</p> <p>ⓘ Press the ESC key if a parameter is to be retained. Correction function for parameters or values (prior to ENTER)</p>	
<p>Level 2</p>	<p>Parameter level</p>	
<p>CONTROL DYNAMIC</p> <p>← →</p> <p>▼ ▲</p>	<p>DYNAMIC APER ✱</p> <p>or</p> <p>DYNAMIC NORM</p>	<ul style="list-style-type: none"> The parameter flashes, switch by pressing ▼ ▲ and ← →
<p>CONTROL xP INT</p> <p>← →</p> <p>▼ ▲</p>	<p>xP INT 16 ✱</p> <p>0.1 ... 99.9</p>	<ul style="list-style-type: none"> The parameter flashes, set by pressing ▼ ▲ + ← →
<p>CONTROL TN INT</p> <p>← →</p> <p>▼ ▲</p>	<p>TN INT 100 ✱</p> <p>3 ... 9999</p>	<ul style="list-style-type: none"> The parameter flashes, set by pressing ▼ ▲ + ← →
<p>CONTROL TV INT</p> <p>← →</p>	<p>TV INT 5 ✱</p> <p>0 ... 999</p>	<ul style="list-style-type: none"> The parameter flashes, set by pressing ▼ ▲ + ← →

7.4.1. Dynamik internal

CONTROL DYNAMIC

This parameter affects the temperature sequence only in case of internal control.

Factory setting: APER (aperiodic)

Possible parameters:

NORM Allows for reaching the setpoint faster – with setpoint change or ramp function – but overshooting of up to 5 % is possible.

APER Ramp function: the increase of temperature occurs temporarily offset and achieves the target temperature without overshooting.
Setpoint change: The temperature increases at the same rate, the target temperature is achieved without overshooting.

ⓘ With both settings constant temperature is achieved after approximately the same time.

7.4.2. Control parameters – XP-, TN-, TV- INTERNAL

In most cases the control parameters preset in the factory are adequate for achieving an optimum temperature sequence.
The control parameters allow adjustment to special control processes..

```
XP INT
  16
```

Setting range: 0.1 ... 99.9

Proportional range >Xp<

The proportional range is the range below the setpoint in which the control circuit reduces the heating capacity from 100% to 0 %

```
TN INT
 100
```

Setting range: 3 ...9999

Reset time >Tn< (Integral component)

Compensation of the remaining control deviation due to proportional regulation. An insufficient reset time may cause instabilities. Excessive reset times will result in unnecessary prolongation of compensation of the control difference.

```
TV INT
  5
```

Setting range: 0 ... 999


Lead time >Tv< (Differential component)

The differential component reduces the transient time. An insufficient lead time will prolong the time required for compensation of disturbance effects and cause high overshooting during run-up. An excessive lead time could cause instabilities (oscillations)

Optimization instructions for the PID control parameters

Optimum setting


Control parameters XP-, TN-, TV- INTERNAL


The heat-up curve reveals possible faulty settings of the control parameter.

Inappropriate settings may produce the following heat-up curves:


Xp too low


Tv/Tn too low

Xp too high or Tv too high

Tv/Tn too high or Xp too high

7.5. MENU SERIAL - BAUDRATE, HANDSHAKE, PARITY

<p>Menu level 1</p> <p>▼ ▲</p> <p>MENU</p> <p>SERIAL</p> <p>↵ ↓</p>	<p>For communication between circulator and a PC or a superordinated process control system the interface parameters of both units must be identical.</p> <p>① For remote control refer to page 97</p> <p>Factory settings: 4800 Baud even hardware handshake</p>	
<p>Level 2</p>	<p>Parameter level</p>	<p>① Press the ESC key if a parameter is to be retained.</p>
<p>SERIAL</p> <p>BAUDRAT ↵ →</p> <p>▼ ▲</p>	<p>BAUDRAT</p> <p>4800 ✱</p> <p>BAUDRAT</p> <p>9600</p> <p>BAUDRAT</p> <p>19200</p> <p>BAUDRAT</p> <p>38400</p>	<ul style="list-style-type: none"> The parameter flashes, switch by pressing ▼ ▲ and ↵
<p>SERIAL</p> <p>PARITY ↵ →</p> <p>▼ ▲</p>	<p>PARITY</p> <p>EVEN ✱</p> <p>PARITY</p> <p>ODD</p> <p>PARITY</p> <p>NO</p>	<ul style="list-style-type: none"> The parameter flashes, switch by pressing ▼ ▲ and ↵ <p>even: Data bits = 7; Stop bits = 1 odd: Data bits = 7; Stop bits = 1 no: Data bits = 8; Stop bits = 1</p>
<p>SERIAL</p> <p>HSHAKE ↵ →</p>	<p>HSHAKE</p> <p>HARD ✱</p> <p>HSHAKE</p> <p>SOFT</p>	<ul style="list-style-type: none"> The parameter flashes, switch by pressing ▼ ▲ and ↵ <p>Xon/Xoff-protocol (Software handshake) Protocol RTS/CTS (Hardware handshake)</p>


7.6. MENU ATC - Absolut Temperature Calibration

MENU
ATC


ATC serves to compensate a temperature difference that might occur between circulator and a defined measuring point in the bath tank because of physical properties.

Example:

1-point calibration

3-point calibration


T_T 1 = Original curve

Principle:

For ATC calibration, in steady state the bath temperature at the location of the temperature sensor (CT) is determined at the respective adjusted working temperature. This value is then set on the circulator in the menu >ATCalibration< under menu item > CTEMP X <.

This can be a 1-point, 2-point or 3-point calibration.

- M = Temperature measuring instrument with temperature sensor
- B = Bath tank (INTernal)
- T = circulator
- CT = Temperature on measuring point
- TT = Temperature on Thermostat

Menu-Ebene 1

ⓘ Press the **ESC** key if parameter is to be retained. Correction function for parameters or values (prior to ENTER).

Level 2	Parameter level	
<p>ATC <I> STATUS</p> <p>Navigation arrows: left, right, up, down</p>	<p>STATUS YES</p> <p>or</p> <p>STATUS NO</p>	<ul style="list-style-type: none"> • The parameter flashes, switch by pressing and ⓘ >NO< Carry out an ATC calibration ⓘ >YES< return to standard operation after calibration.

Menu functions

Level 2	Parameter level	
<pre> ATC <I> TYPE ▼ ▲ </pre>	<pre> TYPE I-POINT ※ TYPE 2-POINT TYPE 3-POINT </pre>	<ul style="list-style-type: none"> The parameter flashes, switch by pressing ▼ ▲ and ← ⓘ A >1-point<, >2-point< or >3-point< calibration can be carried out. <p>The selected calibration is indicated on level 2 by 1 or 2 or 3.</p>
<pre> ATC <I> TMPVAL1 ▼ ▲ ATC <I> CALVAL1 ▼ ▲ </pre>	<pre> TMPVAL1 8000 ESC CALVAL1 79.70 ※ </pre>	<p>The value >TMPVAL< is only indicated</p> <ul style="list-style-type: none"> ⓘ In addition the measured temperature value >CALVAL X< is saved during the next step. Integer digits flash, set by pressing ▼ ▲ + ← Decimal digits flash, set by pressing ▼ ▲ + ←
<p>ⓘ If only a 1-point calibration is carried out, the following menu items are not indicated anymore</p>		
<pre> ATC <I> TMPVAL2 ▼ ▲ </pre>	<pre> TMPVAL2 12000 ESC </pre>	<p>The value is only indicated</p>
<pre> ATC <I> CALVAL2 ▼ ▲ </pre>	<pre> CALVAL2 119.50 ※ </pre>	<ul style="list-style-type: none"> Integer digits flash, set by pressing ▼ ▲ + ← Decimal digits flash, set by pressing ▼ ▲ + ←
<p>ⓘ If only a 2-point calibration is carried out, the following menu items are not indicated anymore</p>		
<pre> ATC <I> TMPVAL3 ▼ ▲ </pre>	<pre> TMPVAL3 16000 ESC </pre>	<p>The value is only indicated</p>
<pre> ATC <I> CALVAL3 ▼ ▲ </pre>	<pre> CALVAL3 159.30 ※ </pre>	<ul style="list-style-type: none"> Integer digits flash, set by pressing ▼ ▲ + ← Decimal digits flash, set by pressing ▼ ▲ + ←

7.6.1. ATC STATUS - YES / NO

```
ATC (I)
STATUS
```

```
STATUS
YES
```

```
STATUS
NO
```

In the second submenu the ATC function for the temperature sensor selected above is activated >YES< or deactivated >NO<.

>YES< (factory setting) The controller of the circulator uses the original curve of the temperature sensor or the new curve measured during the ATC calibration.

Important: Set to >NO< during the calibration process

>NO< An ATC calibration is to be carried out.

Important: Set to >YES< after calibration.

❶ In the > ATC STATUS < >YES< the ATC calibration always affects the current working temperature; also the one set via interface.

7.6.2. CALIBRATION TYPE: 1 -/ 2 -/ 3 POINT

```
ATC (I)
TYPE
```

```
TYPE
1-POINT
```

```
TYPE
2-POINT
```

```
TYPE
3-POINT
```

```
TMPVAL1 CALVAL1
8000 79.70
```

```
TMPVAL2 CALVAL2
12000 119.50
```

```
TMPVAL3 CALVAL3
16000 159.30
```

A >1-point<, >2-point< or >3-point< **calibration** can be carried out.

First geometrically define the location for calibration (measuring point CT), then determine the temperature values of the calibration points.

The type of calibrations also determines the number of the following pairs of values indicated on the VFD-DISPLAY.

Pairs of values:

TMPVAL X: Circulator temperature 1 or 2 or 3 (actual value TT)

The actual temperature of the bath is simultaneously saved with the "calibration value" >CALVAL< and can be indicated for control purposes (value does not flash).

CALVAL X: Calibration temperature 1 or 2 or 3 (actual value CT)

The „calibration value“ is determined with a temperature measuring device and saved under menu item >CALVAL<.

(value flashes )

7.6.3. Example: 3-point calibration for internal control.


In the temperature range from 80 °C to 160 °C the calibration curve of the temperature sensor (TT) is to be adjusted to the actual temperatures at measuring point (CT).

1. Set working temperature setpoint – SETPNT:

Refer to „Direct temperature setting “ page 73

- By pressing the key the circulator switches to the active >SETPOINT< see example on the left: >SETPNT / 1 25.00°C<. The integer digits flash (Example: <25>).
- Change the value to 80.00 °C by pressing the keys and and confirm by pressing the key . The decimal digits flash. Confirm once more by pressing the key .
- The bath is heated up. Wait for approx. 5 minutes until the temperature is constant.

2. Reading of temperature measuring device

Read the value of measuring point CT on the device and enter under menu item >CALVAL X< by using the keypad.

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Menu functions

Level 2	Parameter level	
<pre> ATE <I> TMPVAL1 ←→ ▼ ▲ ATE <I> CALVAL1 ←→ ▼ ▲ </pre>	<pre> TMPVAL1 8000 ESC </pre> <pre> CALVAL1 79.70 ✱ </pre>	<p>The value >TMPVAL< is only indicated In addition the measured value >CALVAL X< is saved during the following step</p> <ul style="list-style-type: none"> Integer digits flash, set by pressing ▼ ▲ (79) + ← Decimal digits flash, set by pressing ▼ ▲ (70) + ← <p>The first of 3 points is calibrated.</p>
Return to 2. Set working temperature value SETPNT: 120.00 °C		
<pre> ATE <I> TMPVAL2 ←→ ▼ ▲ ATE <I> CALVAL2 ←→ ▼ ▲ </pre>	<pre> TMPVAL2 12000 ESC </pre> <pre> CALVAL2 119.50 ✱ </pre>	<p>The value is only indicated</p> <ul style="list-style-type: none"> Integer digits flash, set by pressing ▼ ▲ (119) + ← Decimal digits flash, set by pressing ▼ ▲ (50) + ← <p>The second of 3 points is calibrated.</p>
Return to 2. set working temperature value SETPNT: 160.00 °C		
<pre> ATE <I> TMPVAL3 ←→ ▼ ▲ ATE <I> CALVAL3 ←→ </pre>	<pre> TMPVAL3 16000 ESC </pre> <pre> CALVAL3 159.30 ✱ </pre>	<p>The value is only indicated.</p> <ul style="list-style-type: none"> Integer digits flash, set by pressing ▼ ▲ (159) + ← Decimal digits flash, set by pressing ▼ ▲ (30) + ← <p>The 3-point calibration is completed</p>
<p>4. Return to standard operation</p>		
<pre> ATE <I> STATUS ←→ </pre>	<pre> STATUS YES ✱ OK </pre>	<ul style="list-style-type: none"> Set >YES< after calibration. (Standard operation)

7.7. MENU LIMITS

Menu- Level 1		
		<p>i Press the ESC key if parameter is to be retained. Correction function for parameters or values (prior to ENTER).</p>
Level 2	Parameter level	
		<ul style="list-style-type: none"> Integer digits flash, set by pressing Decimal digits flash, set by pressing
		<ul style="list-style-type: none"> Integer digits flash, set by pressing Decimal digits flash, set by pressing
		<ul style="list-style-type: none"> The value flashes, set by pressing
		<ul style="list-style-type: none"> The value flashes, set by pressing mit

7.7.1. Limits for internal control

SETPOINT MAX / MIN – Maximum and minimum setpoint

Restriction of the adjustable temperature range

The limitation of the operating temperature range effects the temperature setting in the menu with the key **T**.

Only setting of working temperatures which lie within the determined limits is possible

Existing settings for SETPNT 1, -2, -3, as well as those for >OVERTMP< and >SUBTMP < (refer to page 75), are automatically deferred into the limit range.

Setting range: -94,90 °C ... +200,0 °C

SET MAX > SET MIN

Interchange of values is not possible.

Set maximum heating / cooling

The heating and cooling capacity of the unit are adjustable. 100 % corresponds to the technical specification of the equipment.

Setting range:

HEAT MAX – 0 to 100 % in 1 % steps

COOLING MAX – 0 to 100 % in 1 % steps

8. Troubleshooting guide / error messages

-OFF-
ALARM
CODE 01

 
+
WARNING
CODE 03

Alarm with complete shutdown:

If one of the following failures occur a complete, all-pole shutdown of the heater and circulating pump is effected.

„“ lights up and a continuous signal sounds.

The code for the cause of alarm is indicated on the VFD COMFORT-DISPLAY.

Alarm without shutdown:

The code for the cause of alarm is indicated on the VFD COMFORT-DISPLAY. The warning signal sounds in regular intervals. The messages appear every 10 seconds.

Press the key  to stop the signal

ALARM
CODE 01

Error message with ticker: >LOW LEVEL ALARM-FILL MEDIUM <

Low level alarm

The circulator is operated without or insufficient bath fluid.

Switch the unit off with the mains switch, refill bath fluid and switch on!

Tube breakage has occurred (insufficient filling level of bath fluid caused by pumping-out)

Replace the tubing and refill bath liquid.

The float is defect (e.g. transport damage).

Repair by authorized JULABO service personnel.

ALARM
CODE 02

Error message with ticker:

> REFRIGERATOR ALARM-CHECK CONNECTION <

During the self-test after switch-on a short –circuit is detected between pin 2 and pin 4 of the control line or the control line was disconnected during operation.

Reconnect the control line or repair short-circuit.

WARNING
CODE 03
ALARM
CODE 03

Error message with ticker:

> EXCESS TEMPERATURE WARNING-CHECK LIMITS <

Excess temperature warning

or

Excess temperature alarm

Type of warning: set to >warning< or >alarm< (refer to page 75)

WARNING
CODE 04
ALARM
CODE 04

Error message with ticker:

> LOW TEMPERATURE WARNING-CHECK LIMITS <

Low temperature warning

or

Low temperature alarm.

Typ of warning: set to >warning< or >alarm< (refer to page 75)

ALARM
CODE 05

Error message with ticker:

> WORKING SENSOR ALARM-CALL SERVICE <

Cable of working temperature sensor is disconnected or short-circuited.

ALARM
CODE 06

Error message with ticker:

>SENSOR DIFFERENCE ALARM-CHECK VISCOSITY AND PUMP STAGE<

Defect of working or excess temperature protector.

Working temperature sensor and excess temperature protector report a temperature difference of more than 35°C.

ALARM
CODE 07

Error message with ticker:
> INTERNAL HARDWARE ERROR-CALL SERVICE <
Other errors

ALARM
CODE 12


Error in A/D converter

ALARM
CODE 14

Error message with ticker:
> EXCESS TEMPERATURE PROTECTOR ALARM-CHECK ADJUSTMENT <
Excess temperature protector defect.
The protection temperature is below the set working temperature setpoint.
Set the protection temperature to a higher value.

ALARM
CODE 33

Error message with ticker:
> SAFETY SENSOR ALARM-CALL SERVICE <
The cable of the excess temperature protector has been disconnected or short-circuited

By quickly switching off and restarting the unit the alarm is cancelled.
If the error occurs once more after the restart, a remote diagnosis is required.


C-Err
PRESS
OK

Error message with ticker:
> CONFIGURATION ERROR-PRESS OK<
The configuration of the circulator does not correspond with its current application.
Press the  key for a non-recurring, automatic change of the configuration.

In this case please call the Technical Service or an authorized dealer.

Disturbances that are not indicated.

The electronic pump motor is overload-protected by an electronic current limiter. If viscosity of the bath fluid is or becomes too high, the motor stops running.

Mains circuit breakers (resettable) 15 A


9. Electrical connections

Notice

Use shielded cables only.

The shield of the connecting cable is electrically connected to the plug housing.

RS232 serial interface

This port can be used to connect a computer with an RS232 cable for remote control of the circulator.

Beschaltung:

Pin 2	RxD	Receive Data
Pin 3	TxD	Transmit Data
Pin 5	0 V	Signal GND
Pin 7	RTS	Request to send
Pin 8	CTS	Clear to send

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
Connection of flow-through cooler CK300 / CK310

The flow-through cooler **CK300 / CK310** is connected to pins 1 and 3. For the connection use the shielded control cable supplied with the cooler.


Circuit:	Operation	= relay powered
	Alarm	= relay not powered

Pin assignment:

Pin	Signal
1	+24 V (I max. current 25 mA)
2	0 V
3	Alarm relay
4	Reserved - do not use!
5	Cooling pulse

10. Remote control

10.1. Setup for remote control

- Check the interface parameters for both interfaces (on circulator and PC) and make sure they match.
(Serial interface refer to page 87)
- In the menu > MENU CONFIG < set the menu item > REMOTE < to > ON < .
- Connect both units with an interface cable..

Like all parameters which can be entered through the keypad, interface parameters are stored in memory even after the thermostat is turned off.

10.2. Communication with a PC or a superordinated data system

If the thermostat is put into remote control mode via the configuration level, the VFD COMFORT-DISPLAY will read "R -OFF-" = REMOTE STOP. The circulator is now operated via the computer.

In general, the computer (master) sends commands to the circulator (slave).

The thermostat sends data (including error messages) only when the computer sends a query.

In remote control mode, the start command and all values to be set must be resent by the PC via the interface in case of a power interruption.

AUTOSTART is not possible.

A transfer sequence consists of:

- command
- space (↔; Hex: 20)
- parameter (decimal separation with a period)
- end of file (↵; Hex: 0D)

The commands are divided into **in** and **out** commands.

in commands: retrieve parameters

out commands: set parameters

Important times for a command transmission:

To ensure a safe data transfer, the time gap between two commands should be at least 250 ms.

The circulator automatically responds to an **in** command with a data string followed by a LF (Line Feed). The next command should only be sent after 10 ms.

The out commands are valid only in remote control mode.

Command to set the working temperature > SETPNT 1< to 55.5 °C
out_sp_00 ⇔ 55.5↓

Command to retrieve the working temperature > SETPNT 1< in_sp_00↓

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in commands: Asking for parameters or temperature values to be displayed.

Command	Parameter	Response of circulator
version	none	Number of software version (V X.xx)
status	none	Status message, error message (see page 100)
in_pv_00	none	Actual bath temperature.
in_pv_01	none	Heating power being used (%).
in_pv_03	none	Temperature value registered by the safety sensor.
in_pv_04	none	Setpoint temperature of the excess temperature protection
in_sp_01	none	Working temperature „SETPNT 1“
in_sp_02	none	High temperature warning limit „OverTemp“
in_sp_03	none	Low temperature warning limit „SubTemp“
in_par_01	none	Te - Time constant of the external bath.
in_par_02	none	Si - Internal slope
in_par_03	none	Ti - Time constant of the internal bath.
in_par_05	none	Factor pk/ph0: Ratio of max. cooling capacity versus max. heating capacity
in_par_06	none	Xp control parameter of the internal controller.
in_par_07	none	Tn control parameter of the internal controller.
in_par_08	none	Tv control parameter of the internal controller.
in_mode_01	none	Selected setpoint: 0 = SETPNT 1 1 = SETPNT 2 2 = SETPNT 3
in_mode_05	none	Circulator in Stop/Start condition: 0 = Stop 1 = Start
in_mode_08	none	Adjusted control dynamics 0 = aperiodic 1 = standard

10.4. Status messages

Status messages	Description
00 MANUAL STOP	Thermostat in „OFF“ state.
01 MANUAL START	Thermostat in keypad control mode.
02 REMOTE STOP	Thermostat in „r OFF“ state.
04 REMOTE START	Thermostat in remote control mode.

10.5. Error messages

Error messages	Description
-01 TEMP / LEVEL ALARM	Low liquid level alarm.
-02 REFRIGERATOR ALARM	Control cable of the flow-through cooler is short-circuited or interrupted.
-03 EXCESS TEMPERATURE WARNING	High temperature warning.
-04 LOW TEMPERATURE WARNING	Low temperature warning.
-05 TEMPERATURE MEASUREMENT ALARM	Working temperature sensor short-circuited or interrupted.
-06 SENSOR DIFFERENCE ALARM	Sensor difference alarm. Working temperature and safety sensors report a temperature difference of more than 35 °C.
-07 I²C-BUS ERROR	Internal error when reading or writing the I ² C bus.
-08 INVALID COMMAND	Invalid command.
-09 COMMAND NOT ALLOWED IN CURRENT OPERATING MODE	Invalid command in current operating mode.
-10 VALUE TOO SMALL	Entered value too small.
-11 VALUE TOO LARGE	Entered value too large.
-12 WARNING : VALUE EXCEEDS TEMPERATURE LIMITS	Error in A/D converter.
-14 EXCESS TEMPERATURE PROTECTOR ALARM	Excess temperature protector alarm
-30 CONFIGURATION ERROR: CONFIRM BY PRESSING <OK> ON CIRCULATOR	The configuration of the circulator does not conform to its present use. Press  to automatically perform a single modification of the configuration.
-33 SAFETY SENSOR ALARM	Excess temperature sensor short-circuited or interrupted.

11. Cleaning / repairing the unit

Caution

- Always turn off the unit and disconnect the mains cable from the power source before cleaning the unit.
- Prevent humidity from entering into the circulator.
- Electrical connections and any other work must be performed by qualified personnel only.

Cleaning:

For cleaning the bath tank and the immersed parts of the circulator, use low surface tension water (e.g., soap suds).

Clean the outside of the unit using a wet cloth and low surface tension water.

The Thermostat is designed for continuous operation under normal conditions.

Periodic maintenance is not required.

The tank should be filled only with a bath fluid recommended by Xylem Analytics Germany. To avoid contamination, it is essential to change the bath fluid from time to time.

Repairs:

Before asking for a service technician or returning a circulator for repair, please contact our service station (see backside of this manual).

When returning a unit, take care of careful and adequate packing. Xylem Analytics Germany is not responsible for damages that might occur from insufficient packing.

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Comment contacter Xylem Analytics Germany

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SI Analytics is a trademark of Xylem Inc. or one of its subsidiaries.
Xylem Analytics Germany Sales GmbH & Co.KG
© 2022 Xylem, Inc. Version 221123 M