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
Tel.: +49(0)6131.66.5119

Fax: +49(0)6131.66.5001

Service-Sensors.si-analytics@xylem.com

www.XylemAnalytics.com

Operating instructions for redox combination electrodes from SI Analytics®

Putting into operation

The redox combination electrodes are delivered ready for use. If an irrigation cap is located over the membrane and the diaphragm, please remove it. It contains potassium chloride (KCl) solution 3 mol/l. A KCl solution shortage in the electrolyte chamber of the reference system should be topped up. Topping up is unnecessary in the case of low-maintenance electrodes. With these electrodes, irrigation in electrolyte solution is particularly important.

Checking and measuring

For checking and measuring, the gate of the refilling hole has to be open. To do so, please move the slide gate which is integrated in the electrode head into position 1.

The diaphragm has to immerse into the measurement solution. Low-maintenance redox combination electrodes have a gel electrolyte, therefore they are not equipped with such a closure slide. To check the unit, please immerse the redox combination electrode in the testing solution, for instance L 430 from SI Analytics®, and compare the reading with the indicated value. Please do not pour back used checking solutions! Use fresh checking solutions only!

Maintenance

- Soil on metal sensor and diaphragms will lead to measurement deviations. Calcareous deposits may be removed using diluted mineral acids (for instance diluted hydrochloric acid), for organic contamination please use a suitable solvent. Fats are best removed using tenside solutions, and proteins should be dissolved using hydrochloric pepsin solution (e.g. L 510)

- After measuring and cleaning, please rinse the redox combination electrodes with distilled water. Do not rub them dry, and please do not use them as stirrers.

- To refill and store the redox combination electrode, KCl 3 mol/l (for instance, L 300) is to be used.

- When putting on the watering cap of the liquid electrolyte electrodes, please open the slide feed first and close it afterwards!

- Extreme pH values and temperatures will shorten the life time of the redox combination electrodes and restrict the accuracy of the measurements.

For further information, please refer to our website or our specialist literature.

Xylem Analytics Germany GmbH

Am Achalaich 11
82362 Weilheim
Germany

Tel.: +49(0)6131.66.5119

Fax: +49(0)6131.66.5001

Service-Sensors.si-analytics@xylem.com

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