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Cu 1100 A / PLH Copper Halfcell
CU 60 Copper Combination Electrode

**Accuracy when
going to press**

The use of advanced technology and the high quality standard of our products are the result of a continuous development. This may result in differences between this operating manual and your electrode. Also, we cannot guarantee that there are absolutely no errors in this manual. Therefore, we are sure you will understand that we cannot accept any legal claims resulting from the data, figures or descriptions.

**Warranty
declaration**

In case of proper handling, SI ANALYTICS guarantees the quality of the electrode for one year.

Copyright

© Mainz 2013, SI Analytics GmbH
Reprinting - even as excerpts - is only allowed with the explicit written authorization of SI Analytics GmbH, Mainz.
Printed in Germany.

View Halfcell Cu 1100 A / PLH

View Combination electrode CU 60

Commissioning

Combination Electrode CU 60

1	Remove the protection cap.
2	Pull the closing ring downwards so that the refilling opening of the reference electrode is free.
3	Fill ELY/BR/503 electrolyte solution into the refilling opening.
4	Wipe the shaft using a clean paper towel.
5	Place the combination electrode into deionized water for approx. 10 minutes.
6	Before use, place the combination electrode into diluted standard solution for approx. 5 minutes. In order to achieve a short response time, the concentration of the standard solution should be slightly below the measuring range expected.

The combination electrode is now ready to measure. Always leave the refilling opening of the reference electrode open while measuring.

Halfcell Cu 1100 A / PLH

For measurements using the Cu 1100 A / PLH copper electrode, an additional reference electrode is required. The two electrodes together form a combination electrode.

For connecting the Cu 1100 PLH electrode to the meter you would need the cable L1A, L1BNC or L1R, which has to be adjusted to the connection of the meter.

1	Remove the protection cap.
2	Place the electrode into deionized water for approx. 10 minutes.
3	Before use, place the electrode into diluted standard solution for approx. 5 minutes. In order to achieve a short response time, the concentration of the standard solution should be slightly below the measuring range expected.
4	Fill the R 503 reference electrode with ELY/BR/503 electrolyte solution.

The double electrode is now ready to measure.

Note

The reference electrode must contain enough electrolyte solution, so that the level of the electrolyte solution is clearly higher than that of the test sample during the measurement.

Note

We exclusively recommend using the copper electrodes to determine the end point for titration procedures.

Calibration

Calibrate according to the data given in the operating manual of the measuring instrument and the analysis specification.

Sample preparation

Add 2 % ISA/Ca solution. This sample conditioning solution creates optimum conditions for measuring. It provides a constant ionic strength and similar diffusion potentials at the reference electrode in standard and test sample.

Note If you would like to have more detailed information concerning sample preparation and measuring procedures, SI Analytics provides a large number of application reports for various applications.

Response times

The response time depends on the concentration range. It is:

- several seconds with high concentrations,
- several minutes near the detection limit.

The measured value is stable if the value does not change by more than 0.1 mV within 30 seconds.

Interferences

- Complexes with acetate, citrate, ammonia, amino acids, EDTA etc.
- Mercury and silver ions poison the membrane.
- The membrane is photosensitive concerning the setting of the potential.
- Interfering ions: no disturbances with the following conditions:

$(\text{Cu}^{2+}) \times (\text{Cl}^-)^2$	$(\text{Cu}^{2+}) \times (\text{Br}^-)^2$	(Fe^{3+})
$(1,6 \times 10^{-6} \text{ mol/L NO}_3^-)$	$< 1,3 \times 10^{-12} \text{ mol}^3/\text{L}^3$	$< 0,1 (\text{Cu}^{2+})$

Aging

Please note that every electrode undergoes a natural aging process. The response time increases and the slope decreases with the age of the electrode. The following factors shorten the lifetime considerably:

- Incorrect storing
- Special measuring conditions (e.g. organic solutions, frequent measuring with high concentrations of interfering ions)
- High temperatures
- High changes in temperature

The warranty does not cover damage caused by measuring conditions and mechanical damage.

Maintenance

- For combination electrodes: replenish used reference electrolyte solution.
- In case of increased response times of the electrode perform the following activities:
 - Polish the membrane surface using an SF/K polishing strip. To do so, moisten a polishing strip with deionized water, and with circular movements polish the membrane surface for approx. 30 seconds.
 - Subsequently place the electrode into diluted standard solution for approx. 5 minutes.

Storage


Combination electrode CU 60

- Between two measurements, place the combination electrode into diluted standard solution.
- Overnight up to one week: Push the closing ring upwards to close the refilling opening. Then rinse the combination electrode with deionized water, blot it dry using a clean paper towel and store it in a dry place.
- For more than a week: Remove the electrolyte solution from the reference electrode and rinse the combination electrode with deionized water, blot it dry using a clean paper towel and put on the protection cap. Store the combination electrode in a dry place.

Electrode Cu 1100 A / PLH

- Between two measurements, place the electrode into diluted standard solution.
- For storing overnight or longer, rinse the electrode with deionized water, blot it dry using a clean paper towel and put on the protection cap (closing cap of the glass tube). Store the electrode in a dry place.

Typical calibration line of a copper electrode

Technical data
SI Analytics ISE – Halfcell
Cu 1100 A / PLH/Combination Electrode CU 60

Measuring range $6 \times 10^{-4} \dots 6,350 \text{ mg/L Cu}^{2+}$ ($10^{-8} \dots 0,1 \text{ mol/L Cu}^{2+}$)

Reproducibility $\pm 4 \%$

pH range 2 ... 6 (see section INTERFERENCES)

Temperature range 0 ... 80 °C / up to 100 °C for a short time

Membrane resistance $< 1 \text{ M}\Omega$

Length Cu 1100 A: 170 mm (including 50 mm connection head)
Cu 1100 PLH: 153 mm (including 33 mm connection head)
CU 60: 153 mm (including 33 mm connection head)

Diameter Shaft: 12 mm
Connectionhead: 16 mm

Notizen / Notes

SI Analytics

a xylem brand

SI Analytics GmbH

Hattenbergstr. 10

Tel. +49.(0)6131.66.5111

Fax. +49.(0)6131.66.5001

55122 Mainz

Deutschland, Germany, Allemagne, Alemania

E-Mail: si-analytics@xylem.com

www.si-analytics.com

SI Analytics is a trademark of Xylem Inc. or one of its subsidiaries.

© 2013 Xylem, Inc. Version 130828 M **T814849**