

# The Basics of Photometric Measurement

## Part 1: Principles, Optics, AQA, Test Kits I

# Introduction

**Phos** (Greek) for light => Photometry is a measurement method to analyse (aqueous) solutions by means of a light source.

Light (physical) is a spectra of electromagnetic waves, divided into different ranges: Visible light (white light) ranges from approx. 380 – 780 nm



**WTW photometric range  
(190-1100 nm)**

## Photometric / Colorimetric Analysis:

Determination of substances by their specific colour reaction and light absorbance in dependence of their chemical properties at a specific wavelength.

# Introduction – Light Sources and Optics


Specific Wavelengths are obtained by

## Different light sources

- LEDs ( $\lambda_x$ ) = lowest power consumption, lower light intensity
- Tungsten (white light halogene lamp) for VIS range
- Xenon (UV-VIS) => Flash lamp with long life span
- Deuterium (UV) => special lamp, expensive


## Different optical techniques, such as

- Monochromators
- Polychromators
- Filters
- LED



# Optics: Filter and LED Photometer


Filter photometer with reference beam: photoLab<sup>®</sup> S6/S12



LED<sub>λ</sub> + optical Filter – single beam: photoFlex<sup>®</sup> Series



# Monochromator of photoLab<sup>®</sup> 7000 Series




# Measurement modes

What type of measurement is performed in with a photometer?

=> 3 Measurement modes of photometric analysis and their relation

- 1) **Transmission T(%)**: Ratio of light intensity after cuvette ( $I$ ) and before ( $I_0$ )
- 2) **Absorbance**:  
$$\text{Abs}_\lambda = -\log_{10} (T_\lambda)$$
or „extinction of light“ passing the cuvette
- 3) **Concentration**: quantitative analysis of a substance (mg/l, ppm,...) at a defined wavelength based on a calibration curve

# Transmission measurement



Transmission is the ratio of passed light  $I$  / initial light  $I_0$ :


Transmission measurement is also being used to measure **turbidity at 180° angle** (unit FAU, e.g. for quality control) and for turbidity correction in concentration measurement.

# Absorbance : Concentration Measurement

**Absorbance** = „Extinction of light“:


Each substance has a **specific spectra** with absorbance peak(s).

=> Run spectrum to define maximum or optimal peak = wavelength definition for concentration measurement



**Concentration** measurement:

Measurement at **specific wavelength**, obtained by either matching LED, optical filters from white light or monochromator






# Relation %T : Absorbance : Concentration

**Transmission measurement:**

The transmission of a sample varies **exponentially** with thickness and concentration



**Absorbance measurement:**

Absorbance of a sample is **proportional** to thickness of the sample and concentration

Transmission (T%)

100    10    1    0,1

=> Logarithmic correlation

Absorbance  $A = -\log_{10}(T)$

0    1    2    3

} => linear correlation


Concentration (mg/l)

0    4    8    12

# Relation %T : Absorbance : Concentration

## Lambert-Beer's law

Experiments by BOUGUER (1698–1758) and LAMBERT (1728–1777) showed that the **absorbance is dependent on the thickness** of the absorbing layer of the cell used. The relationship between the absorbance **and the concentration** of the analyte in question was discovered by BEER (1825–1863). The combination of these two natural laws led to the derivation of **Lambert-Beer's law**, which can be described in the form of the following equation:



$$A = \epsilon_{\lambda} \cdot c \cdot d$$

$\epsilon_{\lambda}$  = molar absorptivity, in l/mol x cm

$d$  = Path length of the cell, in cm

$c$  = Concentration of the analyte, in mol/l

Source: Operating Instructions of photoLab® S12, Part 1: General Information ([www.wtw.com](http://www.wtw.com))

# Concentration Measurement


The correlation of absorbance/concentration is determined by setting up a **characteristic calibration curve** for each substance (parameter). The chemical reaction must be known:

**Dilution series with defined concentrations, measured at defined  $\lambda$  and cuvette size (pathlength)**

⇒ Characteristic (calibration) curve

⇒ **Unknown sample concentration** can be „read“ from the curve!

Methods/ programs in photometers contain all data and compute result automatically, including various cuvette sizes. Barcoded test kits additionally call up the respective method=program.



# Method data / Program for each parameter

Programmed data for comfortable concentration measurement are consisting in:

$\lambda$  matching the absorbance for determination

Reagent blank  $E_0$  = coloration of reagent


Slope of calibration curve for calculation

Citation & unit (e.g.  $\text{NO}_3\text{-N}$  mg/l)

Factors for citation & unit switch (e.g.  $\text{NO}_3$ ; mmol/l)


Sample blank (e.g. coloration of the sample) is **not** included!

„ $E_{\text{sample}}$ “ comes on top of  $E_0$ : individually for sample, for small volumes mostly negligible



# Prerequisites for concentration measurement

- Coloured solution contains dissolved dye
- Absorption of light leads to a coloration (complementary to  $\lambda$ )
- Color intensity is correlating with concentration
- The chemical reaction of analyte leads to building or disappearance (e.g. COD 4-40 mg/l) of dye in defined reaction time
- Reaction must be selective for the analyte – no cross reaction with other **disturbing** substances
- The developed dye must be stable for time of measurement => e.g. reading within 10 minutes after reaction time w/o color deterioration (see analytical instructions)



# Instrument Check – AQA

- **Self calibration/check &** warm up time (especially for kinetics and spectra)
- **AutoCheck:** photoLab<sup>®</sup> levels meter vs air in the background
- **Zeroing/Baseline:** correcting the meter to „E<sub>0</sub>“ => especially after transport or in changing conditions (temperature...) meters E<sub>0</sub> „drift“ some mAbs => readings become inaccurate (mostly too high)

- **AQA Tools:**

- Optical or liquid filters
- Color solutions, e.g. PhotoCheck<sup>®</sup>
- Selected unscratched zero cuvettes
- Control standards of the substance



# FAQ Commercial test kits in brief

## Measurement range: (MR)

The range is meter (optics) dependent, the reaction has detection limits. MR-values are reaching approx.  $\pm 2 - 2,5$  Abs (test dependent!)

At the lower end **detection limits and tolerances** of the procedure have the biggest **influence on accuracy** of readings: Limitation of chemical procedure, confidence interval and accuracy are often at the lower limit.

=> Scratches, pipette faults etc. **additionally** affect the accuracy of readings! Readings at the lower end **become more inaccurate!**

**Measure in the middle of the MR, if possible!**

Test: A6/25 Ammonium (WTW)

### 2. Measuring range and numt

Measuring range
0.20 - 8.00 mg/l NH <sub>4</sub> -N
0.26 - 10.30 mg/l NH <sub>4</sub> <sup>+</sup>

a): Low end 0,20 mg/l  
confidence interval:  $\pm 0,20$  mg/l

Characteristic data of the procedure:

Sensitivity: Absorbance 0.010 A corresponding (mg/l NH <sub>4</sub> -N)	0.04
Accuracy of a measurement value (mg/l NH <sub>4</sub> -N)	max. $\pm 0.20$

b): Accuracy:  $\pm 0,20$  mg/l

# FAQ – The importance of “photometric Zero”

Performing a zero (see manuals!)

**LED meters, e.g. pHotoFlex® Series**

Portable meters require a zero due to changing conditions, transportation and optics.

**Filter photometers, e.g. photoLab® S12**

In lab with stable conditions, slow drift and often stabilizing reference beam requires less zeroing.

**Spectrophotometers, e.g. photoLab® 7000 Series**

Zero/base line is required for many functions of spectral tasks, concentration mode is similar to filter photometers with reference beam

Influence of meter drift: => Zeroing

Package leaflet of many tests show **sensitivity** by correlating absorbance A (E=A) to mg/l.

Influence can be seen directly:

For COD test 14560, 4-40 mg/l COD, an absorbance of 10 mE means 0.4 mg/l COD. **10 mE drift** without zeroing means 0.4 mg/l or **10% evitable miss-reading** in the low end!

Characteristic data of the procedure:

	1.14560.	1.14540.
Number of lots	23	27
Sensitivity 0.010 E (absorbance) $\hat{=}$ mg/l COD	0.4	2
Accuracy of a measurement value (mg/l COD)	max. $\pm$ 1.8	max. $\pm$ 6