

Determination of Acid Number by Color-Indicator Titration (ASTM D 974)

Description

The determination described in this report is based on the standard ASTM D 974.

Instruments

| | |
|----------------------|-------------------------|
| Titration | TL 7000 or higher |
| Interchangeable unit | WA 10 |
| Electrode | OptiLine 6 |
| Stirrer | Magnetic stirrer TM 235 |
| Lab accessory | Glass beaker 150 ml |
| | Magnetic stirring bars |
| | Balance |

Reagents

| | |
|--|--------------------------------------|
| 1 | KOH in 2-propanol, 0.1 mol/L |
| 2 | Toluene |
| 3 | 2-Propanol |
| 4 | Deionized water |
| 5 | p-Naphtholbenzein indicator solution |
| 6 | Potassium acid phthalate |
| 8 | Phenolphthalein indicator solution |
| All reagents should be of reagent grade or better. | |


Standardization of Potassium Hydroxide Solution

The standardization can be carried out either potentiometric or colorimetric.

If a pH electrode is available. Please refer to the standard method Titer KOH and the application note "Titer strong bases".

The colorimetric indication is possible with phenolphthalein indicator solution. Weigh to the nearest 0.1 mg approximately 0.1 g of potassium acid phthalate, which has been dried for at least 1 h at 110 °C and dissolve in 80-100 mL of water, free of CO₂. Add six drops of the phenolphthalein indicator solution and titrate with linear titration steps to the inflection point. For highest accuracy repeat the titer determination three times and calculate the mean value using the statistic function of the titrator

See the titration curve and method below.



| | | | |
|-------------------------|---------------------|----------------------|-----------|
| Default method | - | | |
| Method type | Automatic Titration | | |
| Modus | Linear | | |
| Measured value | mV-E, 520 nm | | |
| Measuring speed / drift | User defined | Minimum holding time | 05 s |
| | | Maximum holding time | 15 s |
| | | Measuring time | 3 s |
| | | Drift | 10 mV/min |
| Initial waiting time | 10 s | | |
| Linear steps | 0.05 ml | | |
| Damping | None | Titration direction | decrease |
| Pretitration | 6.00 ml | Delay time | 30 s |
| End value | Off | | |
| EQ | Off | Slope value | - |
| Max. titration volume | 12 ml | | |
| Dosing speed | 100% | Filling speed | 30 s |

Calculation:

$$T \text{ [mol/l]} = \frac{W * F2}{(EQ - B) * M * F1}$$

| | | |
|-----|--------|--|
| B | 0 | Blank value |
| W | man | Weight of the sample [g] |
| F2 | 1000 | Conversion factor ml - l |
| EQ1 | | Consumption of titrant until first Equivalence point |
| M | 204,22 | Molecular mass |
| F1 | 1 | Conversion factor |

We recommend to store the exact concentration T into the exchangeable Unit (WA) automatically.

Titration solvent


Prepare by mixing toluene, water, and anhydrous isopropyl alcohol in the ratio 100/1/99.

Indicator solution

Prepare a solution of p-naphtholbenzein in titration solvent equal to 10 g/L.

Blank value determination

Perform a blank titration on 100 ml titration solvent and 0.5 mL of the p-Naphtholbenzein indicator solution.



| | | | |
|-------------------------|---------------------|----------------------|-----------|
| Default method | - | | |
| Method type | Automatic Titration | | |
| Modus | Linear | | |
| Measured value | mV-E, 625 nm | | |
| Measuring speed / drift | User defined | Minimum holding time | 10 s |
| | | Maximum holding time | 40 s |
| | | Measuring time | 4 s |
| | | Drift | 20 mV/min |
| Initial waiting time | 10 s | | |
| Linear steps | 0.02 ml | | |
| Damping | None | Titration direction | decrease |
| Pretitration | off ml | Delay time | 0 s |
| End value | Off | | |
| EQ | Off | Slope value | - |
| Max. titration volume | 0.3 ml | | |
| Dosing speed | 100% | Filling speed | 30 s |

Calculation:

$$\text{Result ml} = EQ1$$

| | |
|-----|---|
| EQ1 | Consumption of titrant at the first Equivalence point |
|-----|---|

The result is stored as global variable M01.


Sample titration

The sample is weight in into a 150 ml beaker. The amount of sample weight is adjusted to the expected acid value.

| Acid number (new and Light oils) | Size of sample [g] | Sensitivity of Weighing [g] |
|----------------------------------|--------------------|-----------------------------|
| 0 – 3 | 20.0 +/- 2.0 | 0.05 |
| Over 3 to 25 | 2.0 +/- 0.2 | 0.01 |
| Over 25 to 250 | 0.2 +/- 0.02 | 0.001 |
| Used or Dark-Colored Oils | | |
| 0 to 25 | 2.0 +/- 0.2 | 0.01 |
| Over 25 to 250 | 0.2 +/- 0.02 | 0.001 |

To the sample 100 ml of solvent and 0.5 ml of the indicator are added. Swirl until the sample is completely dissolved. The color of the indicator must be yellow-orange. If the mixture assumes a green or green-black color then the sample does not have any acid number and the base number can only be determined.

Immerse the OptiLine 6 and the titration tip (if not already happened) and start the titration method.



| | | | |
|-------------------------|---------------------|----------------------|-----------|
| Default method | - | | |
| Method type | Automatic Titration | | |
| Modus | Linear | | |
| Measured value | mV-E, 625 nm | | |
| Measuring speed / drift | User defined | Minimum holding time | 10 s |
| | | Maximum holding time | 40 s |
| | | Measuring time | 4 s |
| | | Drift | 20 mV/min |
| Initial waiting time | 10 s | | |
| Linear steps | 0.04 – 0.05 ml | | |
| Damping | None | Titration direction | decrease |
| Pretitration | off ml | Delay time | 0 s |
| End value | Off | | |
| EQ | Off | Slope value | - |
| Max. titration volume | 6 ml | | |
| Dosing speed | 100% | Filling speed | 30 s |

Maybe use the manual stop!

Calculation:

$$\text{Result [mg KOH/g]} = \frac{(EQ1 - B) * T * M * F1}{W * F2}$$

| | | |
|-----|------|---|
| B | M01 | Blank value |
| EQ1 | | Consumption of titrant at the first Equivalence point |
| T | WA | Actual Concentration of the titrant |
| M | 56.1 | Molecular weight |
| W | man | Sample weight in g |
| F1 | 0.1 | Conversion factor 1 |
| F2 | 1 | Conversion factor 1 |

Any questions? Please contact the application team:

Xylem Analytics Germany Sales GmbH & Co. KG, SI Analytics
Hattenbergstraße 10
D-55122 Mainz, Germany
Telefon: + 49 6131 66 5126
Fax: + 49 6131 66 5101
E-Mail: titration@si-analytics.com

SI Analytics
a xylem brand

Xylem Analytics Germany Sales GmbH & Co. KG · Hattenbergstr. 10 · D-55122 Mainz ·
Germany Telefon: +49 6131.66. 5111 · E-Mail: Info.si-analytics@Xyleminc.com · www.si-analytics.com

Alle Namen sind eingetragene Handelsnamen oder Warenzeichen der Xylem Inc. oder eines seiner Tochterunternehmen. Technische Änderungen vorbehalten.
© 2018 Xylem Analytics Germany Sales GmbH & Co. KG.