



NH 500/2

Ammonia gas-sensitive electrode



a xylem brand

**Note**

The latest version of the present operating manual can be found on the Internet under www.WTW.com.

Copyright


© 2015, WTW GmbH

Reprinting - even as excerpts - is only allowed with the explicit written authorization of WTW GmbH.

Printed in Germany.

1 Overview

Structure



- 1 Protective cap for plug-in connector
- 2 Plug-in connector
- 3 Shaft
- 4 Membrane cap (or storing cap on delivery)
- 5 Ag/AgCl reference electrode
- 6 pH glass electrode
- 7 Gas-permeable membrane

2 Safety

Authorized use

The authorized use of the NH 500/2 consists of its use as an ammonia or ammonium electrode in the laboratory and field.

3 Commissioning

Scope of delivery

- 1 electrode, provided with protective cap and storing cap
- 3 membrane caps
- 1 bottle NH₃ electrolyte solution (50 ml)
- Operating manual

Getting the sensor ready for measuring

On delivery, the electrode is equipped with the storing cap (without the white, gas-permeable membrane).

- Remove the storing cap.
- Rinse the electrode with deionized water.
- Fill approx. 1 ml NH_3 electrolyte solution into a membrane cap.
- Remove air bubbles in the electrolyte by knocking.
- Screw the membrane cap on the electrode.
- Connect the electrode to the meter using the connection cable.
- Keep the storing cap. It serves as a protective cap if the electrode is stored for a longer period of time.

4 Measuring / operation


4.1 Calibration



Typical calibration line of an ammonia electrode

Note

For calibration, please refer to the operating manual of the measuring instrument.



4.2 Sample preparation

Add 1 % alkaline reagent (MZ/ NH_3 /CN) to the aqueous sample. It is essential to observe the WTW application reports and analysis specifications respectively.

4.3 Response time

The response time depends very much on the concentration of the measured ion, the condition of the electrode and the direction of the concentration change.

For a 10-fold concentration change from 50 to 5 mg/l $\text{NH}_4\text{-N}$, the response time t_{99} is < 8 min.

4.4 Interferences

- Volatile bases, e.g. amines.
- Coatings on the membrane (e.g. metal hydroxides)

4.5 Storage

Storage duration	Way of storing
Up to approx. two days	Put in NH_3 electrolyte solution with the membrane cap screwed on
For more than two days	Store the electrode as follows: <ul style="list-style-type: none"> ● Unscrew the membrane cap ● Rinse the electrode ● Fill some drops of tapwater in the storing cap and screw it on. Thus the glass membrane will remain conditioned.

4.6 Aging

The glass part of the electrode that is visible when exchanging the membrane cap is a special pH electrode. Please note that it undergoes a natural aging process. If the required slope cannot be achieved despite having changed the membrane cap, the pH electrode is exhausted.

The pH electrode is unsuitable for conventional pH measurement. This case as well as using unsuitable electrolytes and mechanical damage invalidates any warranty claim.

5 Cleaning and maintenance

5.7 Cleaning

Contamination	Cleaning procedures
Water-soluble contamination	Immerse in deionized water for 10 minutes
Metal hydroxides	Immerse in 10% citric acid

After cleaning, thoroughly rinse with deionized water and recalibrate if necessary.

5.8 Exchanging the membrane cap

The membrane cap is a wear part and has to be replaced from time to time, depending on demand and requirements. A reduced slope, an extended response time or a restricted measuring range indicate that a replacement is necessary. The exchange is carried out according to the chapter, COMMISSIONING.

6 Wear parts and accessories

Wear parts and maintenance means

Description	Model	Order no.
Accessory set, comprising: – 3 membrane caps – 50 ml NH ₃ electrolyte solution	ZBK/NH3/2	180 100

Accessories

Description	Model	Order no.
Connection cable for electrodes with plug-in connector (DIN connector)	AS/DIN	108 110
Connection cable for electrodes with plug-in connector (BNC connector)	AS/BNC	108 114
Alkaline reagent, 10 mol/l NaOH, bottle with 250 ml	MZ/NH ₃ /CN	150 130
Standard solution, 10 g/l ammonium (NH ₄ Cl) , bottle with 1 l	ES/NH ₄	120 240




Note

For further accessories, refer to the WTW catalog or the Internet.

7 Technical data

Dimensions (in mm)



Materials

Shaft, storing cap	POM
Plug-in connector	PPS
Membrane cap	POM
Membrane	PP
Materials with sample contact	POM copolymer, NBR (O ring), PP

Measurement conditions

Measuring ranges at 20 °C	$10^{-6} \dots 5 \cdot 10^{-2} \text{ mol/l NH}_4^+$ $0.02 \dots 900 \text{ mg/l NH}_4^+$
Temperature range	0 ... 50 °C
Depth of immersion	min. 5 mm max. 50 mm
Operating position	vertical, inclined to max. 45 °
Max. allowed overpressure	$< 5 \cdot 10^4 \text{ Pa (0.5 bar)}$

Characteristic data on delivery

Reproducibility	$\pm 2 \%$
-----------------	------------

Xylem |'zīləm|

- 1) The tissue in plants that brings water upward from the roots;
- 2) a leading global water technology company.

We're 12,500 people unified in a common purpose: creating innovative solutions to meet our world's water needs. Developing new technologies that will improve the way water is used, conserved, and re-used in the future is central to our work. We move, treat, analyze, and return water to the environment, and we help people use water efficiently, in their homes, buildings, factories and farms. In more than 150 countries, we have strong, long-standing relationships with customers who know us for our powerful combination of leading product brands and applications expertise, backed by a legacy of innovation.

For more information on how Xylem can help you, go to www.xyleminc.com



Wissenschaftlich-Technische Werkstätten GmbH

Dr.-Karl-Slevogt-Straße 1
D-82362 Weilheim
Germany

Tel: +49 881 183-0
Fax: +49 881 183-420
E-Mail: info@wtw.com
Internet: www.wtw.com