**BEDIENUNGSANLEITUNG
OPERATING MANUAL
MODE D'EMPLOI
INSTRUCCIONES DE
OPERACIÓN**

ba75899defs02 12/2017

ADA S7/IDS

a xylem brand

ADA S7/IDS

| | |
|---|-----------|
| Bedienungsanleitung | 3 |
| Allgemeines | 5 |
| Technische Daten | 6 |
| Inbetriebnahme, Messen, Kalibrieren | 7 |
| Lagerung | 7 |
| Operating Manual | 9 |
| General information | 11 |
| Technical data | 12 |
| Commissioning, measuring, calibration | 13 |
| Storage | 13 |
| Mode d'emploi | 15 |
| Généralités | 17 |
| Caractéristiques techniques | 18 |
| Mise en service, mesure, calibration | 19 |
| Stockage | 19 |
| Instrucciones de operación..... | 21 |
| Información general | 23 |
| Especificaciones técnicas | 24 |
| Puesta en funcionamiento, medición, calibración | 25 |
| Almacenamiento | 25 |

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ADA S7/IDS

IDS pH ADAPTER FOR pH ELECTRODES WITH S7 PLUG-IN CONNECTOR

a xylem brand

Note

The latest version of the present operating manual can be found on the Internet under www.WTW.com.

Copyright

© 2017 Xylem Analytics Germany GmbH
Printed in Germany.

General information

The ADA S7/IDS adapter forms an IDS pH sensor in conjunction with a pH electrode with S7 plug-in connector.

Automatic adapter recognition

The measuring electronics with the stored adapter data is in the adapter head of the ADA S7/IDS. The data include, among other things, the adapter type and series number. With each calibration, the calibration data is written in the adapter and the calibration history is recorded. The data is recalled by the meter when the adapter is connected and is used for measurement and for measured value documentation. Storing the calibration data in the adapter ensures that the correct slope and asymmetry are automatically used if the sensor is operated with several meters.

The digital transmission technique guarantees the failure-free communication with the meter even with long connection cables. If the sensor firmware is enhanced by WTW, it can be updated via the meter.

Note

More detailed information on the automatic adapter recognition or sensor recognition is given in the operating manual of the meter.

Calibration data and calibration history

The connected meter uses for measurement the calibration data stored in the adapter, regardless of the mounted electrode.

Note

Therefore, you have to calibrate the IDS pH sensor each time the electrode with the S7 plug-in connector has been exchanged. Otherwise, the IDS pH sensor will work with the calibration data of the previously used electrode. This could lead to incorrect measured values. When evaluating the calibration history, please note that in between calibrations, the electrode might have been exchanged.


Technical data

| | | | | |
|---|------------------------------------|--|---------------------------------|-----------------------------|
| General data | Electrodes that can be connected | All electrodes with S7 plug-in connector | | |
| | Application characteristics | pH measuring range | Depending on the used electrode | |
| | Allowed temperature range | 0 ... 80 °C (32 ... 176 °F) | | |
| | Typical application | Applications with special pH electrodes | | |
| Characteristics of the IDS measuring technique | Measured parameter | Measuring range | Resolution | Accuracy (± 1 digit) |
| | pH | -2.000 ... 20.000 | 0.001 | ± 0.004 |
| | U [mV] | -1000.0 ... +1000.0 | 0.1 | ± 0.2 |
| Material, electrical connection | Material of the plug enclosure | PPE/PS | | |
| | Combination electrode connection | S7 plug | | |
| | Meter connection | Fixed cable with digital plug | | |
| Connection cable | Length | 1.5 m | | |
| | Diameter | 4.3 mm | | |
| | Smallest allowed bend radius | Fixed installation: | 20 mm | |
| | | Flexible use: | 60 mm | |
| | Plug type | Socket, 4 pins | | |
| Electrical data of the electrode connection | Input resistance | > 5 * 10 ¹² ohm | | |
| | Input current | < 1 * 10 ⁻¹² A | | |

Commissioning, measuring, calibration

Commissioning

Connect the S7 plug-in connector of the pH electrode with the S7 plug of the adapter.

Note

The further commissioning steps are given in the operating manual of the pH electrode.

- Connect the electrode to the meter.
- Calibrate the electrode according to the operating manual of the meter and observe the following rules while doing so:

Calibration and measurement

Connect the adapter with the pH electrode to the meter. All further details are given in the operating manuals of the pH electrode and meter.

Caution: Heed the depth of immersion of the pH electrode. The ADA S7/IDS adapter must not be submersed!

Storage

For storage, leave the adapter plugged on the pH electrode if you do not intend to use both devices otherwise. The storing method depends on the pH electrode. More details are given in the operating manual of the pH electrode.

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What can Xylem do for you?

We're a global team unified in a common purpose: creating innovative solutions to meet our world's water needs. Developing new technologies that will improve the way water is used, conserved, and re-used in the future is central to our work. We move, treat, analyze, and return water to the environment, and we help people use water efficiently, in their homes, buildings, factories and farms. In more than 150 countries, we have strong, long-standing relationships with customers who know us for our powerful combination of leading product brands and applications expertise, backed by a legacy of innovation.

For more information on how Xylem can help you, go to xylem.com.

Service address:

Xylem Analytics Germany
Sales GmbH & Co. KG
WTW
Dr.-Karl-Slevogt-Str. 1
82362 Weilheim
Germany

Tel.: +49 881 183-325
Fax: +49 881 183-414
E-Mail wtw.rma@xylem.com
Internet: www.WTW.com

Xylem Analytics Germany GmbH
Dr.-Karl-Slevogt-Str. 1
82362 Weilheim
Germany

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Copyright © 2017 Xylem Analytics Germany GmbH