**BEDIENUNGSANLEITUNG
OPERATING MANUAL
MODE D'EMPLOI
INSTRUCCIONES DE
OPERACIÓN**

ba75969defs02 08/2015

ADA 94/IDS DIN / ADA 94/IDS BNC

a xylem brand

ADA 94/IDS DIN / ADA 94/IDS BNC

Bedienungsanleitung	2
Allgemeines	2
Inbetriebnahme, Messen, Kalibrieren	3
Technische Daten	4
Operating Manual	5
General information	5
Commissioning, measuring, calibration	6
Technical data	7
Mode d'emploi	8
Généralités	8
Mise en service, mesure, calibration	9
Caractéristiques techniques	10
Instrucciones de operación	11
Generalidades	11
Puesta en funcionamiento, medición, calibración	12
Especificaciones técnicas	13

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General information

The ADA 94/IDS adapter enables the connection of an analog pH electrode (with or without temperature sensor) to the digital meters inoLab Multi 9620, Multi 9630, Multi 9420 and Multi 9430.

Two models of the adapter are available:

Connection of the analog electrode	Adapter
DIN	ADA 94/IDS DIN
BNC	ADA 94/IDS BNC

Automatic adapter recognition

The measuring electronics with the stored adapter data is in the ADA 94/IDS adapter. The data include, among other things, the adapter type and series number. With each calibration, the calibration data is written in the adapter and the calibration history is recorded. The data are called up by the meter and used for measurement and documentation of measured values.

If the adapter firmware is enhanced by WTW, it can be updated via the meter.

More detailed information on the automatic adapter recognition or sensor recognition is given in the operating manual of the meter.

Calibration data and calibration history

The meter uses for measurement the calibration data stored in the adapter, regardless of the mounted electrode.



Therefore, you always have to recalibrate after exchanging the pH electrode. Otherwise, the IDS pH sensor will work with the calibration data of the previously calibrated electrode. This could lead to incorrect measured values. When evaluating the calibration history, please note that in between calibrations, the electrode may have been exchanged.

Commissioning, measuring, calibration

Commissioning

- | | |
|---|--|
| 1 | Unscrew the 4 screws of the cover plate (the screw driver is included in the scope of delivery). |
| 2 | Insert the ADA 94/IDS adapter in the open mounting slot. |
| 3 | Attach the adapter to the meter with the 4 screws. |

Example Multi 9630 with mounted adapter ADA 94/IDS DIN

The further commissioning steps are given in the operating manual of the pH electrode.

- | | |
|---|---|
| 4 | Connect the electrode to the socket of the adapter. |
| 5 | Calibrate the electrode according to the operating manual of the meter. |

All further details are given in the operating manuals of the pH electrode and meter.

Technical data

General data	Electrodes that can be connected	All pH electrodes with DIN or BNC plug head connector
---------------------	----------------------------------	---

Application characteristics	pH measuring range	Depending on the electrode used
	Allowed temperature range	0 ... 80 °C (32 ... 176 °F)
	Typical application	Applications with special pH electrodes

Measuring characteristics of the IDS measuring electronics in the ADA 94/IDS DIN ADA 94/IDS BNC adapter	Measured parameter	Measuring range	Resolution	Accuracy (± 1 digit)
	pH	-2.000 ... 20.000	0.001	± 0.004
	U [mV]	-1000.0 ... +1000.0	0.1	± 0.2
	T [°C]		0.1	± 0.2

Material, electrical connection	Material of the plug enclosure	ABS
	Electrode connection	DIN or BNC
	Meter connection	Mounting slot with contacts

Electrical data of the electrode connection	Input resistance	> 5 * 10 ¹² ohm
	Input current	< 1 * 10 ⁻¹² A

The latest version of the present operating manual can be found on the Internet under www.WTW.com.

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What can Xylem do for you?

We're a global team unified in a common purpose: creating innovative solutions to meet our world's water needs. Developing new technologies that will improve the way water is used, conserved, and re-used in the future is central to our work. We move, treat, analyze, and return water to the environment, and we help people use water efficiently, in their homes, buildings, factories and farms. In more than 150 countries, we have strong, long-standing relationships with customers who know us for our powerful combination of leading product brands and applications expertise, backed by a legacy of innovation.

For more information on how Xylem can help you, go to xyleminc.com.

Service address:

Xylem Analytics Germany
Sales GmbH & Co. KG
WTW
Dr.-Karl-Slevogt-Str. 1
82362 Weilheim
Germany

Tel.: +49 881 183-325
Fax: +49 881 183-414
E-Mail wtw.rma@xyleminc.com
Internet: www.WTW.com

Xylem Analytics Germany GmbH
Dr.-Karl-Slevogt-Str. 1
82362 Weilheim
Germany

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Copyright © 2017 Xylem Analytics Germany GmbH